GlassBONE™

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GB – INSTRUCTION SHEET

DESCRIPTION

GlassBONE™ is an absorbable and synthetic bioceramic composed of natural elements (Ca, Na, Si, P, O). GlassBONE™ is intended to fill temporarily bone defects of traumatic, pathologic or surgical origin until bone regeneration. When implanted in living tissue, GlassBONE™ has the property to form on its surface a calcium phosphate layer, which is substantially equivalent in composition and structure to the hydroxyapatite from the mineral phase of human bone.

REF	size	Volume
GB1.3/5	1 – 3 mm	5 cc
GB1.3/10	1 – 3 mm	10 cc
GB1.3/16	1 – 3 mm	16 cc

INDICATIONS

Filling bone defects in orthopaedic and cranio maxillofacial surgery.

Main indications are:

- Filling bone defects (after tumour resection, after removing iliac bone,...)
 - Vertebral Arthrodesis
 - Prosthesis revision
 - Traumatology
 - Tibial Osteotomy
- GlassBONE™ must not be implanted in osseous sites which can let particles go into articular cavities.
- GlassBONE™ must not be implanted in areas which are subjected to mechanical stresses without osteosynthesis material.

CONTRA-INDICATIONS

As for bone graft in general, GlassBONE™ should not be used in patients with:

- Acute or chronic infection at the surgical site.
- Metabolic diseases (diabetes, hyperparathyroidism, osteomalacia...).
- Severe hepatic and renal dysfunction.
- Any existing condition or disease that will interfere with good bone and tissue healing.

ADVERSE REACTIONS

Possible complications are the same as to be expected of autologous bone grafting procedures: post operative infection, allergic reactions due to the components of product, osseous fracture.

PRECAUTIONS BEFORE USE

- Check the use-by date and check that the packaging is intact before use. Devices for which the packaging is open or damaged or for which the use-by date has elapsed should not be used.
- Control each device before use in order to detect any possible defect. In case of defect, the implant should not be used.
- No studies have been conducted on pregnant women or on breastfeeding-women. For safety reasons, the implantation of GlassBONE™ is inadvisable during pregnancy and breastfeeding.
- GlassBONE™ is intended for single use and should never be reused. Reuse is liable to give rise to pathogen transmission risks and impairment of the physicochemical properties of the implant.

INFORMATION FOR THE PATIENT

- The patient should be informed by the surgeon of the potential risks and undesirable effects resulting from the implantation. His agreement is needed for the surgery to be performed.
- The surgeon should inform the patient receiving this device that the success of the implantation depends on his behaviour.
- The patient should inform his surgeon of any event which would compromise the proper integration of the implant and should undergo regular post-operative examinations.

PREPARATION

- GlassBONE™ could be used alone or mixed with normal saline solution, autologous bone, bone-marrow or patient blood.
- In order to facilitate osseous regeneration, GlassBONE™ should be implanted preferentially in a cancellous bone. In order to enhance the host site vascularization, heighten the host bone walls.
- Remove the excess of particles at the implantation site before suturing.
- Combining GlassBONE™ with any medication during surgery is carried out under the surgeon's responsibility.
- The Good Aseptic Practices and the general principles for drug treatment of the patient must be observed while using GlassBONE™.
- Standard postoperative practices for the treatment and rehabilitation associated with bone grafting must be strictly followed.

STERILISATION AND PACKAGING

- Sterile device packed in a double-b blister and sterilized by ethylene oxide. Sterility is guaranteed until the indicated expiration date, if the package is intact and not opened.
- Re-sterilisation is not allowed.

This product is to be handled or implanted by trained qualified persons having read these instructions for use.

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