REF	
ACT-MS1520	15 x 20 mm
ACT-MS1525	15 x 25 mm
ACT-MM2030	20 x 30 mm
ACT-ML3040	30 x 40 mm

NORAKER

INNOVATIVE BIOMATERIALS

	<div>FR Fabricant</div> <div>EN Manufacturer</div> <div>IT Produttore</div> <div>DE Hersteller</div> <div>ES Fabricante</div> <div>PT Fabricante</div>
---	--

FR

Activioss™

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EN

Activioss™

Resorbable synthetic dental membrane (PLGA) for guided tissue regeneration in periodontal or dental implant surgery.

DESCRIPTION

Activioss™ is a membrane composed of bio-resorbable polymers synthesised from polylactic and polyglycolic acid obtained using a controlled and standardised manufacturing process.

It has a two-layer structure: the fibrous surface, facing the bone, promotes the growth of cells responsible for osteogenesis due to its natural collagen-mimicking architecture. The denser, smooth surface facing the soft tissues prevents the latter from penetrating inside the bone defect.

PROPERTIES

The fibrous layer of Activioss™ absorbs fluids, and blood in particular. Mechanically, this property promotes its stability and adhesion to the walls of the bone defect and, biologically it encourages bone regeneration.

The dense layer has a smooth, shiny appearance. It provides the membrane with a certain amount of rigidity. This enables it to retain its shape, even in the presence of humidity, and makes it easier to handle for prolonged periods in indications in which it is combined with bone substitute biomaterials. Experimental data show that the membrane fulfils its barrier function for a period of over four weeks, before being completely resorbed by the body within four to six months.

The synthetic origin of the polymer, along with its biodegradation into lactic and glycolic acids naturally present in the body mean that the Activioss™ membrane presents a very high level of biocompatibility, demonstrated by experimental data obtained in accordance with the standards currently in force.

INDICATIONS

Activioss™ is a resorbable synthetic membrane designed for use in periodontal and dental implant surgery in the following indications:

> Preservation of alveolar bone volume following extraction

> Covering of bone defects in the event of immediate or delayed implant placement.

> Alveolar crest reconstruction

The Activioss™ membrane may be used alone or combined with a graft or bone filling substitute.

CONTRAINDICATIONS

Activioss™ membrane must not be used in the event of active tissue infection at the implant site.

INSTRUCTIONS FOR USE

Opening the packaging: Activioss™ membrane is packed in a double packaging. Peel open the first pouch and place the second pouch on a sterile surgical drape. Carefully remove the membrane from the pouch using sterile surgical tweezers. Gently speaking, the Activioss™ membrane must be handled with care.

The general principles for sterile handling and good surgical practices must be followed when applying the Activioss™ membrane:

- a full-thickness mucoperiosteal flap and wide release incisions at the level of the teeth bordering the site to be covered must be performed.

- the bone surface under the membrane must be thoroughly cleaned to remove all soft tissue debris.

- the Activioss™ membrane is cut with scissors to adapt its size to the dimensions of the defect exposed. The membrane should overlap the walls of the bone defect by at least two millimetres to fulfil its barrier function.

- position the fibrous layer (matt) facing the bone and the dense layer (shiny) facing the soft tissues. Do not separate the two layers of the membrane.

- slip the membrane under the periosteum at the most fixed part of the flap, then apply Activioss™ over the bone defect. Absorption of blood exudates by the fibrous layer promotes adhesion to the bone surface. At this stage, a sterile spatula may be used to optimise the shape or position of the membrane.

- the high level of mechanical resistance of the Activioss™ membrane allows it to be fixed in place; this may be indicated to prevent it moving under the effect of a load or mobilisation. It may be fixed in place with sutures or pins if required, but this is not essential.

-while the membrane is exposed, limit contact with saliva as much as possible in order to avoid any bacterial contamination.

- the mucoperiosteal flap should be sutured, totally covering the membrane (for example using simple sutures and deep mattress sutures).

- if synthetic or biological materials are used to fill the bone defect, it should not be excessively filled.

POSTOPERATIVE CARE

- The patient should receive regular follow-up and be prescribed a mouthwash for two weeks following the surgical procedure. Level I analgesics are recommended for forty-eight hours in the event of pain. Oral and dental hygiene must be maintained. It is recommended that the site be brushed using a post-surgical toothbrush for forty-eight hours following the surgery. Healthcare workers will inform the patient about the necessity to keep away from heavy smoking and to have proper oral and dental hygiene for two weeks following the surgery.

- If the membrane is exposed, the dehiscence generally heals of its own accord in a few weeks.

- Removal of the membrane is not generally necessary, except in the case of acute infection. If removal of the membrane is necessary, the adjacent tissues must be anaesthetized using a local anaesthetic and an incision must be made on the part immediately adjacent to the residual membrane. After carefully folding back the tissue surrounding the remaining portion, the membrane can be excised and the curettage performed to remove inflamed or infected tissue.

- To allow good bone regeneration under the membrane, do not perform any surgical re-entry for 4 to 6 months following the surgery.

PRECAUTIONS FOR USE

Activioss™ membrane should only be used by dentists and surgeons trained in its use.

To ensure the success of any surgical treatment inside the mouth, it is necessary to eliminate any underlying bacterial infection and ensure good oral hygiene. Consequently, before any surgery is performed, patients must receive hygiene treatment based on instructions related to oral hygiene and periodontal cleaning, including dental scaling or root planning and occlusal adjustment, if indicated.

If endo-osseous implants are used, only use the membrane in combination with a stable implant and not with a view to obtaining the primary stability of the implant.

The Activioss™ membrane must not be used in multiple layers.

The Activioss™ membrane must not be used as a bone filling material.

The Activioss™ membrane must not be imbedd with active substances (antibiotics, growth factors, bisphosphonates). The Activioss™ membrane must not be used on the skin and integuments.

Particular caution is required when using Activioss™ membrane in the following patients:

- patients with poorly controlled metabolic disorders (e.g. diabetes, osteomalacia, thyroid disorders);
- patients receiving concomitant treatments, such as long-term corticosteroid therapy, anti-rejection drugs, radiotherapy or chemotherapy;
- patients with an auto-immune disease;
- heavy smokers;
- patients with chronic infections (for example osteomyelitis) at the surgical site.

The effects on pregnant or breastfeeding women are not known.

POSSIBLE COMPLICATIONS

The adverse effects observed with other synthetic membranes, such as local inflammation, may occur.

The potential complications that may arise following any surgical procedures may develop.

Oedema at the surgical site, partial necrosis of the flap, haemorrhage, dehiscence, haematoma, increased sensitivity and pain, bone loss, erythema and local inflammation, infection and complications related to the local anaesthetic may potentially develop.

STORAGE AND HANDLING INSTRUCTIONS

The Activioss™ membrane is a sterile product designed exclusively for single use during the course of the same surgical procedure.

Do not re-sterilise Activioss™ membrane.

Do not use the product if the sterile packaging is damaged or open.

Do not use if the membrane is soiled before application.

Do not cut using a heat-generating method.

Do not use after the expiry date.

Store at room temperature (preferably < 25°C) away fr om heat and moisture.

The Activioss™ membrane does not interfere with, and is not affected by a magnetic environment.

The Activioss™ membrane can be disposed of or discarded together with domestic garbage unless it has been in contact with the patient's biological fluids, in which case it should be thrown with biological waste.

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