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GlassBone® Granules

Substitut osseux bioactif

Indications

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Substitut osseux bioactif

Indications

Glassbone Granules is a synthetic and biocompatible bone substitute device (bioactive glass 45S5), intended for the filling, reconstruction and/or fusion of bone defects. It is indicated in adults and children in the event of loss or lack of bone substance, for bone defects of traumatic, pathological or surgical origin, when autologous solutions are not applicable or sufficient, in orthopaedic surgery, cranio-maxillofacial (including ENT) surgery, spinal surgery and dental/oral surgery. It must be used by qualified surgeons (orthopaedists, neurosurgeons, dental surgeons, maxillofacial surgeons, stomatologists) trained in bone grafting and fixation techniques who have read these instructions for use.

Contraindications

- In irradiated bone (according to radiological criteria indicating osteonecrosis)
- In case of acute or chronic infection of the surgical site
- In patients with severe trauma with open wounds susceptible to infection near the defect.
- In patients with known allergies to bioactive glass or its components (Ca²⁺, PO₄³⁻ Na⁺ and Si (OH)₄).
- In patients with pre-existing conditions or diseases that may interfere with good tissue healing (e.g., patients treated with biphosphonates).
- In patients who have undergone or will undergo chemotherapy or radiotherapy at or near the implant site.
- To replace structures subject to high mechanical stress
- In case of chronic metabolic condition (uncontrolled diabetes, hyperparathyroidism, osteomalacia, chronic inflammatory rheumatism, severe or fracturing osteoporosis, etc.).
- In severe renal and hepatic infections.
- In patients receiving treatment known to affect the bones.
- In the neonatal department, in infants weighing less than 3 kg
- In the event of a meningeal breach in spinal/cranial surgery
- No studies have been conducted to date in pregnant women and there are no data on use during breastfeeding. As a security measure, implantation of Glassbone Granules is not recommended during pregnancy and lactation. To date, we have no specific data on the use of Glassbone Granules in dental surgery in the paediatric population (only in oral surgery).

Composition of Glassbone Granules

Glassbone Granules (45S5) consist only of elements found naturally in bone tissue: 24.5% calcium oxide (CaO), 6% phosphate pentoxide (P₂O₅), 24.5% sodium oxide (Na₂O) and 45% silicon dioxide (SiO₂).

Mechanism of action / Performance

During its implantation, Glassbone Granules is in contact with bone and biological fluids. The release of ions during resorption will allow the surface formation of a layer of carbonated hydroxyapatite, whose composition and structure are similar to the mineral phase of bone. This layer gives Glassbone Granules its osteoconduction property and creates a strong chemical link between the granules and living tissue.

Instructions for use

- Check the expiry date. Do not use the product if it is exceeded.
- Check each device before use to detect any deterioration. If deterioration is present, do not use the implant.
- Eliminate all soft and/or pathological tissue from the implantation site.
- Open the outer blister (sterile barrier) and remove the inner blister on the sterile field.
- Once the surgical site has been prepared, open the inner blister.
- If necessary, mix the substitute in a sterile cup with another constituent (physiological saline, autologous bone and/or blood)
- Using a sterile instrument, fill in the defect. Do not compress the material in the site or blot the blood/moisture in the positioned graft.
- After placement of Glassbone Granules, ensure the primary closure of the soft tissues at the graft site. Resorbable or non-resorbable membranes can also be used for the closure. In dental surgery, Glassbone Granules can also be used with a fibrin-rich plasma (FRP) membrane. As with any particle grafting procedure, primary closure of the surgical site is essential.

Warnings and precautions for use

Regarding the surgical procedure

- The general principles of aseptis and patient medication must be observed when using Glassbone Granules.
- The combination of any drug substance with Glassbone Granules during implantation is the responsibility of the surgeon.
- Manipulate Glassbone Granules with a surgical instrument to avoid piercing surgical gloves.
- It is advisable to revive the recipient site before implantation.
- Completely fill the defect with Glassbone Granules. In dental surgery, excess filling is recommended.
- Avoid placing granules outside of the bone defect. Remove them if necessary.
- If it moves/migrates, the bioactive glass can cause wear of the joints and interfere with movement. Prevention of movement and granule migration is essential for proper bone formation.
- Do not apply excessive pressure to the defect. Excessive pressure may cause embolization of fat in the bloodstream.
- Glassbone Granules does not have sufficient mechanical strength to withstand load bearing before hard tissue is formed. When used in load bearing areas such as mandible fractures, standard internal or external stabilization techniques should be followed to achieve rigid stabilization in all planes.
- In dental use, Glassbone Granules should not be used in applications with immediate load bearing. Placing a dental implant under load can be done 5 to 6 months after filling. Immediate loading is the responsibility of the surgeon.
- It is necessary to follow the usual post-operative treatment and rehabilitation procedures associated with bone grafts.
- The closure of the operative site depends on the surgery performed and the surgical site (membrane, sutures, etc.).
- The use of a membrane is recommended in cases where the defect is significant or if there is limited bone retention. Primary closure of defects, preferably without tension, is strongly recommended.

Regarding the medical device

- Glassbone Granules is a device that resorbs over time to make way for regenerated bone. There is currently no clinical study available that demonstrates complete resorption of the granules.
- Glassbone Granules has not been evaluated for safety and compatibility in the MRI environment.
- This device does not harden like cement.
- Glassbone Granules is a sterile disposable device and must never be re-sterilized or reused. Reuse may cause contamination and impairment of bone substitute performance.

Side effects

No side effect directly linked to the device has been reported to date. However, an unknown allergy to one of the constituents of the product, a delay in consolidation or a fusion failure may be possible. Possible complications are the general complications due to surgery or anaesthesia: post-surgical symptoms (pain, redness, inflammation, oedema, haematomas, seroma, swelling, etc.), post-operative infection, delay in consolidation, loss of reduction of the fracture, fusion failure, loss of bone graft, protrusion of the graft. The following complications can also occur in dental use: tooth sensitivity, gingival recession, subsidence of the flap, abscess formation, resorption or anklyosis of the treated root, etc. These complications are the same as those that can occur with autologous bone grafting. Any serious incident which may occur in connection with Glassbone Granules must be notified to NORAKER and to the competent authority of the Member State in which the patient or the surgeon is established.

Informations for the patient

- The patient must be informed by their surgeon of the potential risks and adverse effects of implantation and agree to the proposed procedure.
- The surgeon should inform the patient who is the recipient of this device that the success of implantation depends on their behaviour and good compliance with post-operative hygiene instructions.
- The patient must report any incident to their surgeon that could compromise the proper integration of the implant and undergo postoperative checks.

Sterilization

Glassbone Granules is a single-use device, sterilized with ethylene oxide. Sterility is guaranteed until the expiration date, if the sterile barrier has not been opened or damaged. The allowable limit of residual EtO defined by the manufacturer is 0.012 mg/MD.

Storage & Disposal

The devices must be stored in their original unopened packaging in a clean, dry place, away from direct sunlight and at a recommended temperature between 15°C and 25°C. Disposal of the device should be carried out in accordance with local regulations and practices.

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