# GlassBone®

## Injectable Putty

FR – NOTICE D'UTILISATION - Substitut osseux bioactif .....	6
EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute.....	8
IT – ISTRUZIONI PER L'USO - Sostituto osseo bioattivo .....	10
DE – GEBRAUCHSANWEISUNG - Bioaktives Knochenersatzmaterial .....	12
ES - INSTRUCCIONES DE USO - Sustituto óseo bioactivo .....	15
PT - INSTRUÇÕES DE UTILIZAÇÃO - Substituto de osso bioativo .....	17
BG – ИНСТРУКЦИИ ЗА УПОТРЕБА - Биоактивен заместител на костна тъкан .....	19
CS – NÁVOD K POUŽITÍ - Bioaktivní kostní náhrada .....	22
EL – ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Βιοενεργό υποκατάστατο οστών .....	24
HR – UPUTE ZA UPORABU - Bioaktivni koštani nadomjestak .....	26
NL – GEBRUIKSAANWIJZING - Bioactieve botvervanger .....	28
PL – INSTRUKCJA OBSŁUGI - Bioaktywny substytut kości .....	31
SL – NAVODILA ZA UPORABO - Bioaktivni nadomestek kosti .....	33
SR – УПУТСТВО ЗА УПОТРЕБУ - Биоактивна замена за кости .....	35
-AR - تعليمات الاستخدام - بديل عظام نشط حيويًا -	37
LV – LIETOŠANAS INSTRUKCIJA – Bioaktīvs kaulu aizstājējs .....	39

REF	VOL
GB-IP1.0	1cc
GB-IP2.5	2.5cc
GB-IP5.0	5cc
GB-IP10	10cc

**NORAKER®**  
THE BIOGLASS® COMPANY


### NORAKER

60, Avenue Rockefeller

69008 LYON, France

Tel : +33 (0)4 78 93 30 92, Fax : +33 (0)4 72 35 94 37

[contact@noraker.com](mailto:contact@noraker.com), [www.noraker.com](http://www.noraker.com)

	FR	Conserver au sec
	EN	Keep dry
	IT	Tenere all'asciutto
	DE	Vor Feuchtigkeit schützen
	ES	Manténgase seco
	PT	Manter seco
	BG	Да се пази сухо
	CS	Uchovávejte v suchu
	EL	Να διατηρείται στεγνό
	HR	Čuvati na suhom
NL	Droog houden	
PL	Trzymać w suchym miejscu	
SL	Hranite na suhem	
SR	Држати на сувом	
AR	يجب الحفاظ على المنتج جافاً	
LV	Glabāt sausumā	
	FR	Conserver à l'abri de la lumière du soleil
	EN	Store away from sunlight
	IT	Tenere lontano dalla luce del sole
	DE	Vor Sonnenlicht schützen
	ES	Mantener alejado de la luz solar
	PT	Manter afastado da luz solar
	BG	Да се пази от слънчева светлина
	CS	Uchovávejte mimo dosah slunečního záření
	EL	Να φυλάσσεται μακριά από το ηλιακό φως
	HR	Držati podalje od sunčeve svjetlosti
NL	Uit de buurt van zonlicht houden	
PL	Trzymać z dala światła słonecznego	
SL	Hranite izven dosega sončne svetlobe	
SR	Држати даље од сунчеве светлости	
AR	يجب الحفاظ على المنتج بعيداً عن أشعة الشمس	
LV	Glabāt prom no saules gaismas	
	FR	Limite de température 15-25°C (stockage)
	EN	Temperature limit 15-25°C (storage)
	IT	Limite di temperatura 15-25 °C (conservazione)
	DE	Temperaturbegrenzung 15-25 °C (Lagerung)
	ES	Límite de temperatura 15-25 °C (almacenamiento)
	PT	Limite de temperatura 15-25°C (armazenamento)
	BG	Температурна граница 15-25°C (съхранение)
	CS	Teplotní omezení 15–25 °C (skladování)
	EL	Όριο θερμοκρασίας 15-25 °C (αποθήκευση)
	HR	Ograničenje temperature: 15-25 °C (skladištenje)
NL	Temperatuurslimiet 15-25°C (opslag)	
PL	Limit temperatury 15–25°C (przechowywanie)	
SL	Temperaturna meja 15-25 °C (shranjevanje)	
SR	Температурна граница 15-25°C (чување)	
AR	حد درجة الحرارة 15-25 درجة مئوية (التخزين)	
LV	Temperatūras ierobežojums 15-25°C (uzglabāšana)	
	AR	Ne pas réutiliser
	EN	Do not reuse
	IT	Non riutilizzare
	DE	Nicht erneut verwenden
	ES	No reutilizar
	PT	Não reutilizar
	BG	Да не се употребява повторно
	CS	Nepoužívejte opakovaně
	EL	Μην επαναχρησιμοποιείτε
	HR	Nemojte ponovno upotrebljavati
NL	Niet hergebruiken	
PL	Nie używać ponownie	
SL	Ponovna uporaba ni dovoljena	
SR	Не користити поново	
AR	لا تقوم بإعادة استخدام المنتج	
LV	Nelietot atkārtoti	
	FR	Stérile, Irradiation gamma
	EN	Sterile, Gamma irradiation
	IT	Sterile, irradiazione gamma
	DE	Sterilisiert durch Gamma-Bestrahlung
	ES	Estéril, radiación Gamma
	PT	Estéril, Irradiação gama
	BG	Стерилно, Гама облъчване
	CS	Sterilní, gama záření
	EL	Αποστειρωμένο με ακτινοβολία γάμμα
	HR	Sterilizirano gama-zračenjem
NL	Steriel, gammastraling	
PL	Produkt sterylny, napromienianie gamma	
SL	Sterilno, gama obsevanje	
SR	Стерилно, гама зрачење	
AR	مُعقم بالتعرض لأشعة جاما	
LV	Sterilizēts ar gamma starojumu	
	FR	Ne pas restériliser
	EN	Do not reesterilize
	IT	Non risterilizzare
	DE	Nicht erneut sterilisieren
	ES	No reesterilizar
PT	Não reesterilizar	

	BG	Да не се стерилизира повторно
	CS	Nesterilizujte opakovaně
	EL	Μην επαναποστειρώνετε
	HR	Nemojte ponovno sterilizirati
	NL	Niet opnieuw steriliseren
	PL	Nie sterylizować ponownie
	SL	Ne sterilizirajte ponovno
	SR	Ne sterilisati ponovo
	AR	لا تقم بإعادة تعقيم المنتج
	LV	Nesterilizēt atkārtoti
	FR	Ne pas utiliser si l'emballage est endommagé et consulter la notice d'utilisation
	EN	Do not use if package is damaged and consult instruction for use
	IT	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	DE	Nicht verwenden, wenn die Verpackung beschädigt ist und Gebrauchsanweisung beachten
	ES	No utilizar si el envase está dañado. Consúltense las instrucciones de uso.
	PT	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	BG	Да не се употребява, ако опаковката е повредена и да се следват инструкциите за употреба
	CS	Pokud je obal poškozen, nepoužívejte a přečtěte si návod k použití
	EL	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευθείτε τις οδηγίες χρήσης
	HR	Nemojte upotrebljavati ako je pakiranje oštećeno; pogledajte upute za uporabu
	NL	Niet gebruiken als de verpakking beschadigd is. Raadpleeg de gebruiksaanwijzing
	PL	Nie używać w razie uszkodzenia opakowania i zapoznać się z instrukcją obsługi
	SL	Ne uporabljajte, če je ovojnina poškodovana. Oglejte si navodila za uporabo
	SR	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu
	AR	لا تستخدم المنتج في حالة تلف العبوة واطلع على تعليمات الاستخدام
	LV	Neļietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju
	FR	Attention Notice d'utilisation
	EN	Instruction for Use
	IT	Attenzione Istruzioni per l'uso
	DE	Vorsicht
	ES	Precaución, Instrucciones de uso
	PT	Instruções de cuidado de utilização
	BG	Внимание Инstrukция за употреба
	CS	Upozornění: Návod k použití
	EL	Προσοχή, συμβουλευτείτε τις οδηγίες χρήσης
	HR	Oprez
	NL	Opgelet gebruiksaanwijzing
	PL	Przeostroga – instrukcja obsługi
	SL	Navodila za uporabo
	SR	Oppez Упутство за употребу
	AR	تنبيه تعليمات الاستخدام
	LV	Lietošanas instrukcija
	FR	Numéro de lot
	EN	Batch number
	IT	Numero di lotto
	DE	Chargenbezeichnung
	ES	Código de lote
	PT	Numero do lote
	BG	Номер на партида
	CS	Číslo dávky
	EL	Αριθμός παρτίδας
	HR	Broj serije
	NL	Partijnummer
	PL	Numer partii
	SL	Številka serije
	SR	Број серије
	AR	رقم المجموعة
	LV	Partijas numurs
	FR	Date limite d'utilisation
	EN	Use by date
	IT	Data di scadenza
	DE	Verwendbar bis
	ES	Fecha de caducidad
	PT	Data de validade
	BG	Срок на годност
	CS	Použitelné do
	EL	Ημερομηνία λήξης
	HR	Rok uporabe
	NL	Houdbaarheidsdatum
	PL	Termin przydatności
	SL	Datum uporabe
	SR	Употребити до датума
	AR	تاريخ انتهاء الصلاحية
	LV	Izlietot līdz (datums)
	FR	Référence
	EN	Reference

REF

IT	Riferimento
DE	Bestellnummer/Artikelnummer
ES	Referencia
PT	Referência
BG	Референция
CS	Reference
EL	Αναφορά
HR	Ref. oznaka
NL	Referentie
PL	Odnosnik
SL	Reference
SR	Референца
AR	المرجع

 0459  
2017

LV	Atsauce
FR	Dispositif marqué CE par l'organisme notifié GMED. Mis sur le marché en 2017
EN	Device marked CE by the GMED notified body. Placed on the market in 2017
IT	Dispositivo marcato CE dall'organismo notificato GMED. Immeso sul mercato nel 2017
DE	Produkt mit CE-Kennzeichnung durch die benannte Stelle GMED. Markteinführung 2017
ES	Producto con marcado CE del organismo notificado GMED. Comercializado desde 2017
PT	Dispositivo com a marca CE emitida pelo organismo notificado pela GMED. Lançado no mercado em 2017
BG	Изделие с ЕС-знак от GMED нотифициран орган. Пуснато на пазара от 2017 г.
CS	Zdravotnický prostředek s označením EC notifikovaným orgánem GMED. Uvedeno na trh v roce 2017
EL	Προϊόν με σήμανση EC από τον κοινοποιημένο οργανισμό GMED. Διατέθηκε στην αγορά το 2017
HR	Prijavljeno tijelo GMED označilo je proizvod oznakom EC. Stavljeno na tržište 2017.
NL	Apparaat EG-gemarkeerd door de GMED-aangemelde instantie. Op de markt gebracht in 2017
PL	Wyrób oznaczony symbolem WE przez jednostkę notyfikowaną GMED. Rok wprowadzenia na rynek: 2017
SL	Priprimoček ES – označen s strani priglašenega organa GMED. Na trg je bil dan leta 2017
SR	Уређај са ознаком ЕС од стране GMED пријављеног тела. Стављен на тржиште 2017
AR	المعتددة. طرح في السوق في عام GMED الجهاز مميز بعلامة المطابقة الأوروبية من قبل جهة 2017
LV	Ierīcei ir CE marķējums, ko veikusi tās pilnvarotā iestāde GMED. Laista tirgū 2017. gadā

FR	Fabricant
EN	Manufacturer
IT	Produttore
DE	Hersteller
ES	Fabricante
PT	Fabricante
BG	Производител
CS	Výrobce
EL	Κατασκευαστής
HR	Proizvođač
NL	Fabrikant
PL	Producent
SL	Proizvajalec
SR	Произвођач
AR	الشركة المصنعة
LV	Ražotājs

NORAKER  
60, Avenue Rockefeller  
69008 Lyon, FRANCE  
Tel : +33 (0)4 78 93 30 92  
Fax : +33 (0)4 72 35 94 37

FR	Date de fabrication
EN	Manufacturing date
IT	Data di produzione
DE	Herstellungsdatum
ES	Fecha de fabricación
PT	Data de fabrico
BG	Дата на производство
CS	Datum výroby
EL	Ημερομηνία κατασκευής
HR	Datum proizvodnje
NL	Productiedatum
PL	Data produkcji
SL	Datum proizvodnje
SR	Датум производње
AR	تاريخ التصنيع

MD

LV	Ražošanas datums
FR	Dispositif Medical
EN	Medical Device
IT	Dispositivo medico
DE	Medizinprodukt
ES	Producto sanitario
PT	Dispositivo Médico
BG	Медицинско изделие
CS	Zdravotnický prostředek
EL	Ιατροτεχνολογικό προϊόν
HR	Medicinski proizvod

	NL	Medisch apparaat
	PL	Wyrób medyczny
	SL	Medicinski pripomoček
	SR	Медицинско средство
	AR	جهاز طبي
	LV	Medicīnas ierīce
	FR	Consulter les instructions d'utilisation
	EN	Consult instructions for use
	IT	Consultare le istruzioni per l'uso
	DE	Siehe Gebrauchsanweisung
	ES	Consultarse las instrucciones de uso
	PT	Consultar as instruções de utilização
	BG	Следвайте инструкциите за употреба
	CS	Přečtěte si návod k použití
	EL	Συμβουλευθείτε τις οδηγίες χρήσης
	HR	Obvezno pročitati upute
	NL	Raadpleeg de gebruiksaanwijzing
	PL	Zapoznać się z instrukcją obsługi
	SL	Ravnajte se po navodilih za uporabo
	SR	Погледати упутства за употребу
	AR	اطلع على تعليمات الاستخدام
	LV	Skatīt lietošanas instrukciju
	FR	Volume (cc / cm <sup>3</sup> )
	EN	Volume (cc / cm <sup>3</sup> )
	IT	Volume (cc / cm <sup>3</sup> )
	DE	Volumen (cm <sup>3</sup> )
	ES	Volumen (cm <sup>3</sup> )
	PT	Volume (cc / cm <sup>3</sup> )
	BG	Обем (cc / cm <sup>3</sup> )
	CS	Objem (cc/cm <sup>3</sup> )
	EL	Όγκος (cc/cm <sup>3</sup> )
	HR	Volumen (cc / cm <sup>3</sup> )
	NL	Volume (cc/cm <sup>3</sup> )
	PL	Pojemność (cc / cm <sup>3</sup> )
	SL	Prostornina (cc/cm <sup>3</sup> )
	SR	Запремина (cc / cm <sup>3</sup> )
	AR	الحجم (سنتيمتر مكعب/سم <sup>3</sup> )
	LV	Tilpums (cc / cm <sup>3</sup> )
	FR	Système de barrière stérile unique avec emballage de protection intérieur
	EN	Single sterile barrier system with protective packaging inside
	IT	Sistema a barriera sterile singola con confezione protettiva interna
	DE	Einfach-Steril-Barrieresystem mit inliegender Schutzverpackung
	ES	Sistema protector estéril único con envase de protección interno
	PT	Sistema de barreira estéril única com embalagem protetor no interior
	BG	Система с единична стерилна бариера с вътрешна защитна опаковка
	CS	Jediný sterilní bariérový systém s ochranným obalem uvnitř
	EL	Ενιαίο σύστημα στείρου φραγμού με εσωτερική προστατευτική συσκευασία
	HR	Sustav s jednostrukom sterilnom barijerom sa zaštitnim unutarnjim pakiranjem
	NL	Enkel steriel barrièresysteem met beschermende verpakking aan de binnenkant
	PL	System pojedynczej bariery sterylnej z wewnętrznym opakowaniem ochronnym
	SL	Sistem sterilne pregrade z zaščitno embalažo v notranjosti
	SR	Једноструки стерилни систем баријере са заштитним паковањем унутра
	AR	نظام حاجز فردي معقم بداخله عبوة واقية
	LV	Viena sterila barjersistēma ar aizsargpakojumu iekšpusē
	FR	Identifiant unique de dispositif
	EN	Unique Device Identifier
	IT	Identificativo univoco del dispositivo
	DE	Einmalige Produktkennung
	ES	Identificador único del dispositivo (UDI)
	PT	Identificador Único do Dispositivo
	BG	Уникален идентификатор на изделието
	CS	Jedinečný identifikátor zařízení
	EL	Μοναδικό αναγνωριστικό προϊόντος
	HR	Jedinstvena identifikacija proizvoda
	NL	Unieke apparaat-ID
	PL	Unikalny identyfikator urządzenia
	SL	Enolični identifikator naprave
	SR	Јединствени идентификатор средства
	AR	معرف الجهاز الفريد
	LV	Ierīces unikālais identifikators
	FR	Mise à jour du document
	EN	Last update
	IT	Ultimo aggiornamento
	DE	Letzte Überarbeitung
	ES	Última actualización
	PT	Última atualização
	BG	Последна промяна
	CS	Poslední aktualizace
	EL	Τελευταία ενημέρωση
	HR	Zadnji put revidirano

**VOL**

**UDI**

**2022/05**

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## **EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute**

### **Description and indications**

GlassBone Injectable Putty is a synthetic, resorbable, biocompatible and bioactive substitute device (bioactive glass 45S5), intended for the filling, reconstruction and / or fusion of bone defects or gaps of the skeletal system, in orthopedic surgery, spine, cranio-maxillofacial surgery and ENT.

It is indicated in case of loss or lack of bone substance for bone defects of traumatic, pathological or surgical origin when autologous solutions are not applicable or sufficient in orthopedics, neurosurgery, cranio maxillo facial and otorhinolaryngology surgery in children and adult population:

- Fusion or reconstruction of deformities and degenerative diseases in spine
- Fusion or reconstruction of deformities and degenerative bone pathologies in orthopedic
- Filling and reconstruction of bone defects due to resection of tumors, cyst or infection and in case of prosthetic revision
- Filling after surgical bone defect (donor sites after removal of autograft, trepanation,...)
- Filling after removal of cholesteatoma
- Filling and reconstruction due to maxilla and periodontium pathologies (in adult only).

It must be applied by qualified surgeons (orthopedists, neurosurgeons, cranio-maxillofacial surgeons, stomatologists and ENT surgeons) trained in bone grafting and fixation techniques and who have read these instructions for use.

The radio-opacity of GlassBone Injectable Putty makes it possible to discern bone substitute granules following their implantation. As the granules resorb, the radio-opacity of the bone defect approaches that of the surrounding bone.

GlassBone Injectable Putty biomaterial is a non-hardening paste, ready to use.

### **Contraindications**

GlassBone Injectable Putty should not be used:

- In case of chronic or acute infection not treated with appropriate therapy.
- In patients who have suffered severe trauma with open external wounds near the defect, which are likely to become infected.
- In patients with a known allergy to bioactive glass or its constituents ( $\text{Ca}^{2+}$ ,  $\text{PO}_4^{3-}$ ,  $\text{Na}^+$  and  $\text{Si}(\text{OH})_4$ ), polyethylene glycol and / or glycerol.
- In patients with pre-existing conditions or disease that may interfere with proper tissue healing.
- In patients who have undergone or will undergo chemotherapy or radiation therapy at or near the site of implantation.
- In the irradiated bone (according to radiological criteria indicating osteonecrosis).
- To replace structures subjected to high mechanical stresses.
- During severe renal and hepatic infections.
- In conjunction with a treatment known to affect the skeleton.
- In case of unsutured meningeal breach in cranio-spinal surgery.
- In neonatology service.
- To date, we do not have any studies conducted in pregnant women or data related to use during breastfeeding.

As a safety measure, the implantation of GlassBone Injectable Putty is not recommended during the periods of pregnancy and lactation.

Composition

GlassBone Injectable Putty is an injectable paste composed of 45S5 bioactive glass granules (69%) mixed with a resorbable binder combining polyethylene glycol (12%) and glycerol (19%).

### **Mechanism of action / Performance / Benefits**

During its implantation, GlassBone Injectable Putty is in contact with bone and biological fluids. Bioactive granules consist only of elements naturally present in bone tissue (Ca, P, Na, Si, O). After implantation, the binder is rapidly absorbed, allowing fluid circulation and cellular and vascular colonization of interstitial spaces between bioactive glass granules.

The resorption of bioactive glass will allow the formation on the surface of a layer of carbonate hydroxyapatite, which composition and structure are similar to the mineral phase of bone. This layer gives the granules their osteoconduction property and makes it possible to create a strong chemical link between the granules and the living tissues. In vitro cell culture tests have also shown that the released ions promote the proliferation, differentiation, and activity of the cells responsible for the formation of bone tissue.

The claimed clinical performance is the filling, reconstruction and / or fusion of bone defects allowing the regeneration of the bone.

The associated benefits are the absence/decrease of autologous bone sampling, a decrease in pain and an improvement in quality of life.

### **Instruction for use**

- Check the expiry date. Do not use the product if it is exceeded.
- Do not use the device if the temperature indicator has turned or if the sterile packaging is damaged.
- Check each device before use, in order to detect any deterioration. If this is the case, do not use the implant.
- GlassBone Injectable Putty does not require mixing, humidification or other preparation prior to use. Humidification could dissolve the binder prematurely.
- Remove any soft and/or pathological tissues from the implantation site.
- Open the outer pouch (sterile barrier) and take out the inner pouch on the sterile field.
- Once the surgical site is prepared, open the internal pouch and remove the applicator.
- Unscrew the wing cap. Depending on the size of the defect and at the convenience of the surgeon, the tip can be unscrewed as well.
- Position the applicator at the level of the bone defect and press the plunger to gently fill the defect. Avoid direct contact between the delivery device and the patient's tissues. Use a sterile instrument if necessary. The device can also be placed directly in an intervertebral cage in spinal surgery.
- Fully fill the defect with GlassBone injectable Putty.
- It is also possible to deposit the paste in a sterile cup and then perform the implantation with a sterile instrument.
- The closure of the surgical site depends on the surgery performed and on the surgical site (membrane, stitches ...).
- This device is MR safe.

### **Warnings and precautions for use**

#### ***In relation to the surgical procedure***

- The general principles of asepsis and patient medication should be respected when using GlassBone Injectable Putty. GlassBone Injectable Putty does not substitute antibiotic therapy treatment during infection.
- The combination of any drug substance with GlassBone Injectable Putty during implantation is the responsibility of the surgeon.
- Handle GlassBone Injectable Putty with a surgical instrument to avoid piercing surgical gloves.
- It is advisable to trim up the recipient site before implantation.
- Avoid placing paste outside the bone defect. Remove it if necessary.
- If positioned outside the implantation site, moving or migrating, bioactive glass can cause wear of the joints and interfere with movement.
- Do not exert excessive pressure on the defect. Excessive pressure could cause an embolism of fat or paste in the bloodstream or cause the paste to be extruded beyond the implantation site, damaging the surrounding tissues.
- GlassBone Injectable Putty does not have sufficient mechanical strength to withstand a load before the bone tissue is formed. When used in load-bearing areas such as mandible fractures, standard internal or external stabilization techniques should be used to achieve rigid stabilization in all planes.
- It is necessary to follow the usual post-operative procedures of treatment and rehabilitation associated with bone grafts.

#### ***In relation to the medical device***

- GlassBone Injectable Putty is a device that resorbs over time to make way for a regenerated bone. The binder is reabsorbed in a few days. Regarding the granules, no clinical studies currently available demonstrate complete resorption.
- GlassBone Injectable Putty is a non-hardening, ready-to-use paste.
- GlassBone Injectable Putty is a sterile single-use device and must not be re-sterilized or reused under any circumstances. Reuse can cause contamination and impaired performance of the bone substitute.

### **Adverse effects**

No side effects directly related to the device have been reported to date. However, an unknown allergy to one of the constituents of the product could occur. Delayed union or failed fusion may also occur depending on the patient's metabolism.

Possible complications are general complications due to surgery or anesthesia: post-surgical symptoms (pain, redness, inflammation, edema, hematomas, seroma, swelling ...), postoperative infection, delay in consolidation, loss of fracture reduction, fusion failure, fracture, loss of bone graft, protuberance of the graft. These complications are the same as those that can occur with autologous bone grafts.

Any serious incident that may occur in connection with GlassBone Injectable Putty must be notified to NORAKER and the competent authority of the Member State in which the patient or surgeon is established.

#### **Patient information**

- The patient must be informed by the surgeon of the potential risks and adverse effects related to implantation and give his/her agreement to the proposed intervention.
- The surgeon must inform the patient receiving this device that the success of the implantation also depends on its behavior and good compliance with the post-operative hygiene instructions.
- The patient must report to his surgeon any event that may compromise the proper integration of the implant and undergo postoperative controls.
- After surgery, an implant card with its fascicle is filled out by the medical staff and given to the patient. He will have to keep it for life. Also, it is advisable to scan it when returning home.
- The SSCP (Summary of Safety and Clinical Performance) of the device is available on the manufacturer's website ([www.noraker.com](http://www.noraker.com)) (or on EUDAMED as soon as available).

#### **Sterilization and packaging**

Sterile device packed in sealed double-pouches and sterilized by gamma radiation (sterile barrier ensured by external pouch). Sterility is guaranteed until the expiry date if the sterile barrier has not been opened or damaged. Re-sterilization is prohibited.

#### **Storage and disposal**

Devices must be stored in their original unopened packaging, in a dry, clean place, protected from sunlight and at a recommended temperature between 15°C and 25°C.

The disposal of the device (packaging, syringe, remaining paste) must be carried out in accordance with local regulations and practices, at the risk of exposing users and patients to pathogens and contaminating the waste circuit.

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