BiologicGlass

Injectable Putty

| | |
|---|----|
| FR – NOTICE D'UTILISATION - Substitut osseux bioactif | 6 |
| EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute..... | 8 |
| IT – ISTRUZIONI PER L'USO - Sostituto osseo bioattivo | 10 |
| DE – GEBRAUCHSANWEISUNG - Bioaktives Knochenersatzmaterial | 12 |
| ES - INSTRUCCIONES DE USO - Sustituto óseo bioactivo | 15 |
| PT - INSTRUÇÕES DE UTILIZAÇÃO - Substituto de osso bioativo | 17 |
| BG – ИНСТРУКЦИИ ЗА УПОТРЕБА - Биоактивен заместител на костна тъкан | 19 |
| CS – NÁVOD K POUŽITÍ - Bioaktivní kostní náhrada | 21 |
| EL – ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Βιοενεργό υποκατάστατο οστών | 23 |
| HR – UPUTE ZA UPORABU - Bioaktivni koštani nadomjestak..... | 26 |
| NL – GEBRUIKSAANWIJZING - Bioactieve botvervanger..... | 28 |
| PL – INSTRUKCJA OBSŁUGI - Bioaktywny substytut kości..... | 30 |
| SL – NAVODILA ZA UPORABO - Bioaktivni nadomestek kosti | 33 |
| SR – УПУТСТВО ЗА УПОТРЕБУ - Биоактивна замена за кости | 35 |
| -AR - تعليمات الاستخدام - بديل عظام نشط حيويًا-AR..... | 37 |
| LV – LIETOŠANAS INSTRUKCIJA – Bioaktīvs kaulu aizstājējs | 39 |

| REF | VOL |
|------------|-------|
| XBG-IP1.0 | 1cc |
| XBG -IP2.5 | 2.5cc |
| XBG -IP5.0 | 5cc |
| XBG -IP10 | 10cc |

NORAKER[®]
THE BIOGLASS[®] COMPANY


NORAKER

60, Avenue Rockefeller

69008 LYON, France

Tel : +33 (0)4 78 93 30 92, Fax : +33 (0)4 72 35 94 37

contact@noraker.com, www.noraker.com

| | | |
|---|---|--|
|  | FR | Conserver au sec |
| | EN | Keep dry |
| | IT | Tenere all'asciutto |
| | DE | Vor Feuchtigkeit schützen |
| | ES | Manténgase seco |
| | PT | Manter seco |
| | BG | Да се пази сухо |
| | CS | Uchovávejte v suchu |
| | EL | Να διατηρείται στεγνό |
| | HR | Čuvati na suhom |
| NL | Droog houden | |
| PL | Trzymać w suchym miejscu | |
| SL | Hranite na suhem | |
| SR | Држати на сувом | |
| AR | يجب الحفاظ على المنتج جافاً | |
| LV | Glabāt sausumā | |
|  | FR | Conserver à l'abri de la lumière du soleil |
| | EN | Store away from sunlight |
| | IT | Tenere lontano dalla luce del sole |
| | DE | Vor Sonnenlicht schützen |
| | ES | Mantener alejado de la luz solar |
| | PT | Manter afastado da luz solar |
| | BG | Да се пази от слънчева светлина |
| | CS | Uchovávejte mimo dosah slunečního záření |
| | EL | Να φυλάσσεται μακριά από το ηλιακό φως |
| | HR | Držati podalje od sunčeve svjetlosti |
| NL | Uit de buurt van zonlicht houden | |
| PL | Trzymać z dala światła słonecznego | |
| SL | Hranite izven dosega sončne svetlobe | |
| SR | Држати даље од сунчеве светлости | |
| AR | يجب الحفاظ على المنتج بعيداً عن أشعة الشمس | |
| LV | Glabāt prom no saules gaismas | |
|  | FR | Limite de température 15-25°C (stockage) |
| | EN | Temperature limit 15-25°C (storage) |
| | IT | Limite di temperatura 15-25 °C (conservazione) |
| | DE | Temperaturbegrenzung 15-25 °C (Lagerung) |
| | ES | Límite de temperatura 15-25 °C (almacenamiento) |
| | PT | Limite de temperatura 15-25°C (armazenamento) |
| | BG | Температурна граница 15-25°C (съхранение) |
| | CS | Teplotní omezení 15–25 °C (skladování) |
| | EL | Όριο θερμοκρασίας 15-25 °C (αποθήκευση) |
| | HR | Ograničenje temperature: 15-25 °C (skladištenje) |
| NL | Temperatuurslimiet 15-25°C (opslag) | |
| PL | Limit temperatury 15–25°C (przechowywanie) | |
| SL | Temperaturna meja 15-25 °C (shranjevanje) | |
| SR | Температурна граница 15-25°C (чување) | |
| AR | حد درجة الحرارة 15-25 درجة مئوية (التخزين) | |
| LV | Temperatūras ierobežojums 15-25°C (uzglabāšana) | |
|  | AR | Ne pas réutiliser |
| | EN | Do not reuse |
| | IT | Non riutilizzare |
| | DE | Nicht erneut verwenden |
| | ES | No reutilizar |
| | PT | Não reutilizar |
| | BG | Да не се употребява повторно |
| | CS | Nepoužívejte opakovaně |
| | EL | Μην επαναχρησιμοποιείτε |
| | HR | Nemojte ponovno upotrebljavati |
| NL | Niet hergebruiken | |
| PL | Nie używać ponownie | |
| SL | Ponovna uporaba ni dovoljena | |
| SR | Не користити поново | |
| AR | لا تقوم بإعادة استخدام المنتج | |
| LV | Nelietot atkārtoti | |
|  | FR | Stérile, Irradiation gamma |
| | EN | Sterile, Gamma irradiation |
| | IT | Sterile, irradiazione gamma |
| | DE | Sterilisiert durch Gamma-Bestrahlung |
| | ES | Estéril, radiación Gamma |
| | PT | Estéril, Irradiação gama |
| | BG | Стерилно, Гама облъчване |
| | CS | Sterilní, gama záření |
| | EL | Αποστειρωμένο με ακτινοβολία γάμμα |
| | HR | Sterilizirano gama-zračenjem |
| NL | Steriel, gammastraling | |
| PL | Produkt sterylny, napromienianie gamma | |
| SL | Sterilno, gama obsevanje | |
| SR | Стерилно, гама зрачење | |
| AR | مُعقم بالتعرض لأشعة جاما | |
| LV | Sterilizēts ar gamma starojumu | |
|  | FR | Ne pas restériliser |
| | EN | Do not reesterilize |
| | IT | Non risterilizzare |
| | DE | Nicht erneut sterilisieren |
| | ES | No reesterilizar |
| PT | Não reesterilizar | |

| | | |
|---|----|---|
| | BG | Да не се стерилизира повторно |
| | CS | Nesterilizujte opakovaně |
| | EL | Μην επαναποστειρώνετε |
| | HR | Nemojte ponovno sterilizirati |
| | NL | Niet opnieuw steriliseren |
| | PL | Nie sterylizować ponownie |
| | SL | Ne sterilizirajte ponovno |
| | SR | Ne sterilisati ponovo |
| | AR | لا تقم بإعادة تعقيم المنتج |
| | LV | Nesterilizēt atkārtoti |
| | FR | Ne pas utiliser si l'emballage est endommagé et consulter la notice d'utilisation |
| | EN | Do not use if package is damaged and consult instruction for use |
| | IT | Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso |
| | DE | Nicht verwenden, wenn die Verpackung beschädigt ist und Gebrauchsanweisung beachten |
| | ES | No utilizar si el envase está dañado. Consúltense las instrucciones de uso. |
| | PT | Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização |
|  | BG | Да не се употребява, ако опаковката е повредена и да се следват инструкциите за употреба |
| | CS | Pokud je obal poškozen, nepoužívejte a přečtěte si návod k použití |
| | EL | Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευθείτε τις οδηγίες χρήσης |
| | HR | Nemojte upotrebljavati ako je pakiranje oštećeno; pogledajte upute za uporabu |
| | NL | Niet gebruiken als de verpakking beschadigd is. Raadpleeg de gebruiksaanwijzing |
| | PL | Nie używać w razie uszkodzenia opakowania i zapoznać się z instrukcją obsługi |
| | SL | Ne uporabljajte, če je ovojnina poškodovana. Oglejte si navodila za uporabo |
| | SR | Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu |
| | AR | لا تستخدم المنتج في حالة تلف العبوة واطلع على تعليمات الاستخدام |
| | LV | Neļietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju |
| | FR | Attention Notice d'utilisation |
| | EN | Instruction for Use |
| | IT | Attenzione Istruzioni per l'uso |
| | DE | Vorsicht |
| | ES | Precaución, Instrucciones de uso |
| | PT | Instruções de cuidado de utilização |
|  | BG | Внимание Инstrukция за употреба |
| | CS | Upozornění: Návod k použití |
| | EL | Προσοχή, συμβουλευτείτε τις οδηγίες χρήσης |
| | HR | Oprez |
| | NL | Opgelet gebruiksaanwijzing |
| | PL | Przeostroga – instrukcja obsługi |
| | SL | Navodila za uporabo |
| | SR | Oppez Упутство за употребу |
| | AR | تنبيه تعليمات الاستخدام |
| | LV | Lietošanas instrukcija |
| | FR | Numéro de lot |
| | EN | Batch number |
| | IT | Numero di lotto |
| | DE | Chargenbezeichnung |
| | ES | Código de lote |
| | PT | Numero do lote |
|  | BG | Номер на партида |
| | CS | Číslo dávky |
| | EL | Αριθμός παρτίδας |
| | HR | Broj serije |
| | NL | Partijnummer |
| | PL | Numer partii |
| | SL | Številka serije |
| | SR | Број серије |
| | AR | رقم المجموعة |
| | LV | Partijas numurs |
| | FR | Date limite d'utilisation |
| | EN | Use by date |
| | IT | Data di scadenza |
| | DE | Verwendbar bis |
| | ES | Fecha de caducidad |
| | PT | Data de validade |
|  | BG | Срок на годност |
| | CS | Použitelné do |
| | EL | Ημερομηνία λήξης |
| | HR | Rok uporabe |
| | NL | Houdbaarheidsdatum |
| | PL | Termin przydatności |
| | SL | Datum uporabe |
| | SR | Употребити до датума |
| | AR | تاريخ انتهاء الصلاحية |
| | LV | Izlietot līdz (datums) |
| | FR | Référence |
| | EN | Reference |

REF

| | |
|----|-----------------------------|
| IT | Riferimento |
| DE | Bestellnummer/Artikelnummer |
| ES | Referencia |
| PT | Referência |
| BG | Референция |
| CS | Reference |
| EL | Αναφορά |
| HR | Ref. oznaka |
| NL | Referentie |
| PL | Odnosnik |
| SL | Reference |
| SR | Референца |
| AR | المرجع |

 0459
2017

| | |
|----|--|
| LV | Atsauce |
| FR | Dispositif marqué CE par l'organisme notifié GMED. Mis sur le marché en 2017 |
| EN | Device marked CE by the GMED notified body. Placed on the market in 2017 |
| IT | Dispositivo marcato CE dall'organismo notificato GMED. Imnesso sul mercato nel 2017 |
| DE | Produkt mit CE-Kennzeichnung durch die benannte Stelle GMED. Markteinführung 2017 |
| ES | Producto con marcado CE del organismo notificado GMED. Comercializado desde 2017 |
| PT | Dispositivo com a marca CE emitida pelo organismo notificado pela GMED. Lançado no mercado em 2017 |
| BG | Изделие с ЕС-знак от GMED нотифициран орган. Пуснато на пазара от 2017 г. |
| CS | Zdravotnický prostředek s označením EC notifikovaným orgánem GMED. Uvedeno na trh v roce 2017 |
| EL | Προϊόν με σήμανση EC από τον κοινοποιημένο οργανισμό GMED. Διατέθηκε στην αγορά το 2017 |
| HR | Prijavljeno tijelo GMED označilo je proizvod oznakom EC. Stavljeno na tržište 2017. |
| NL | Apparaat EG-gemarkeerd door de GMED-aangemelde instantie. Op de markt gebracht in 2017 |
| PL | Wyrób oznaczony symbolem WE przez jednostkę notyfikowaną GMED. Rok wprowadzenia na rynek: 2017 |
| SL | Priprimoček ES – označen s strani priglašenega organa GMED. Na trg je bil dan leta 2017 |
| SR | Уређај са ознаком ЕС од стране GMED пријављеног тела. Стављен на тржиште 2017 |
| AR | المعتدة. طرح في السوق في عام GMED الجهاز مميز بعلامة المطابقة الأوروبية من قبل جهة 2017 |
| LV | Ierīcei ir CE marķējums, ko veikusi tās pilnvarotā iestāde GMED. Laista tirgū 2017. gadā |

| | |
|----|----------------|
| FR | Fabricant |
| EN | Manufacturer |
| IT | Produttore |
| DE | Hersteller |
| ES | Fabricante |
| PT | Fabricante |
| BG | Производител |
| CS | Výrobce |
| EL | Κατασκευαστής |
| HR | Proizvođač |
| NL | Fabrikant |
| PL | Producent |
| SL | Proizvajalec |
| SR | Произвођач |
| AR | الشركة المصنعة |
| LV | Ražotājs |

NORAKER
60, Avenue Rockefeller
69008 Lyon, FRANCE
Tel : +33 (0)4 78 93 30 92
Fax : +33 (0)4 72 35 94 37

| | |
|----|-----------------------|
| FR | Date de fabrication |
| EN | Manufacturing date |
| IT | Data di produzione |
| DE | Herstellungsdatum |
| ES | Fecha de fabricación |
| PT | Data de fabrico |
| BG | Дата на производство |
| CS | Datum výroby |
| EL | Ημερομηνία κατασκευής |
| HR | Datum proizvodnje |
| NL | Productiedatum |
| PL | Data produkcji |
| SL | Datum proizvodnje |
| SR | Датум производње |
| AR | تاريخ التصنيع |

MD

| | |
|----|-------------------------|
| LV | Ražošanas datums |
| FR | Dispositif Medical |
| EN | Medical Device |
| IT | Dispositivo medico |
| DE | Medizinprodukt |
| ES | Producto sanitario |
| PT | Dispositivo Médico |
| BG | Медицинско изделие |
| CS | Zdravotnický prostředek |
| EL | Ιατροτεχνολογικό προϊόν |
| HR | Medicinski proizvod |

| | | |
|--|----|--|
| | NL | Medisch apparaat |
| | PL | Wyrób medyczny |
| | SL | Medicinski pripomoček |
| | SR | Медицинско средство |
| | AR | جهاز طبي |
| | LV | Medicīnas ierīce |
| | FR | Consulter les instructions d'utilisation |
| | EN | Consult instructions for use |
| | IT | Consultare le istruzioni per l'uso |
| | DE | Siehe Gebrauchsanweisung |
| | ES | Consultarse las instrucciones de uso |
| | PT | Consultar as instruções de utilização |
| | BG | Следвайте инструкциите за употреба |
| | CS | Přečtěte si návod k použití |
| | EL | Συμβουλευθείτε τις οδηγίες χρήσης |
| | HR | Obvezno pročitati upute |
| | NL | Raadpleeg de gebruiksaanwijzing |
| | PL | Zapoznać się z instrukcją obsługi |
| | SL | Ravnajte se po navodilih za uporabo |
| | SR | Погледати упутства за употребу |
| | AR | اطلع على تعليمات الاستخدام |
| | LV | Skatīt lietošanas instrukciju |
| | FR | Volume (cc / cm ³) |
| | EN | Volume (cc / cm ³) |
| | IT | Volume (cc / cm ³) |
| | DE | Volumen (cm ³) |
| | ES | Volumen (cm ³) |
| | PT | Volume (cc / cm ³) |
| | BG | Обем (cc / cm ³) |
| | CS | Objem (cc/cm ³) |
| | EL | Όγκος (cc/cm ³) |
| | HR | Volumen (cc / cm ³) |
| | NL | Volume (cc/cm ³) |
| | PL | Pojemność (cc / cm ³) |
| | SL | Prostornina (cc/cm ³) |
| | SR | Запремина (cc / cm ³) |
| | AR | الحجم (سنتيمتر مكعب/سم ³) |
| | LV | Tilpums (cc / cm ³) |
| | FR | Système de barrière stérile unique avec emballage de protection intérieur |
| | EN | Single sterile barrier system with protective packaging inside |
| | IT | Sistema a barriera sterile singola con confezione protettiva interna |
| | DE | Einfach-Steril-Barrieresystem mit inliegender Schutzverpackung |
| | ES | Sistema protector estéril único con envase de protección interno |
| | PT | Sistema de barreira estéril única com embalagem protetor no interior |
| | BG | Система с единична стерилна бариера с вътрешна защитна опаковка |
| | CS | Jediný sterilní bariérový systém s ochranným obalem uvnitř |
| | EL | Ενιαίο σύστημα στείρου φραγμού με εσωτερική προστατευτική συσκευασία |
| | HR | Sustav s jednostrukom sterilnom barijerom sa zaštitnim unutarnjim pakiranjem |
| | NL | Enkel steriel barrièresysteem met beschermende verpakking aan de binnenkant |
| | PL | System pojedynczej bariery sterylnej z wewnętrznym opakowaniem ochronnym |
| | SL | Sistem sterilne pregrade z zaščitno embalažo v notranjosti |
| | SR | Једноструки стерилни систем баријере са заштитним паковањем унутра |
| | AR | نظام حاجز فردي معتم بداخله عبوة واقية |
| | LV | Viena sterila barjersistēma ar aizsargpakojumu iekšpusē |
| | FR | Identifiant unique de dispositif |
| | EN | Unique Device Identifier |
| | IT | Identificativo univoco del dispositivo |
| | DE | Einmalige Produktkennung |
| | ES | Identificador único del dispositivo (UDI) |
| | PT | Identificador Único do Dispositivo |
| | BG | Уникален идентификатор на изделието |
| | CS | Jedinečný identifikátor zařízení |
| | EL | Μοναδικό αναγνωριστικό προϊόντος |
| | HR | Jedinstvena identifikacija proizvoda |
| | NL | Unieke apparaat-ID |
| | PL | Unikalny identyfikator urządzenia |
| | SL | Enolični identifikator naprave |
| | SR | Јединствени идентификатор средства |
| | AR | معرف الجهاز الفريد |
| | LV | Ierīces unikālais identifikators |
| | FR | Mise à jour du document |
| | EN | Last update |
| | IT | Ultimo aggiornamento |
| | DE | Letzte Überarbeitung |
| | ES | Última actualización |
| | PT | Última atualização |
| | BG | Последна промяна |
| | CS | Poslední aktualizace |
| | EL | Τελευταία ενημέρωση |
| | HR | Zadnji put revidirano |

VOL

UDI

2023/08

| | |
|----|------------------------|
| NL | Laatste update |
| PL | Ostatnia aktualizacja |
| SL | Zadnja posodobitev |
| SR | Последње ажурирање |
| AR | تاريخ آخر تحديث |
| LV | Pēdējais atjauninājums |

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EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute

Description and indications

BiologicGlass Injectable Putty is a synthetic, resorbable, biocompatible and bioactive substitute device (bioactive glass 45S5), intended for the filling, reconstruction and / or fusion of bone defects or gaps of the skeletal system, in orthopedic surgery, spine, cranio-maxillofacial surgery and ENT.

It is indicated in case of loss or lack of bone substance for bone defects of traumatic, pathological or surgical origin when autologous solutions are not applicable or sufficient in orthopedics, neurosurgery, cranio maxillo facial and otorhinolaryngology surgery in children and adult population:

- Fusion or reconstruction of deformities and degenerative diseases in spine
- Fusion or reconstruction of deformities and degenerative bone pathologies in orthopedic
- Filling and reconstruction of bone defects due to resection of tumors, cyst or infection and in case of prosthetic revision
- Filling after surgical bone defect (donor sites after removal of autograft, trepanation,...)
- Filling after removal of cholesteatoma
- Filling and reconstruction due to maxilla and periodontium pathologies (in adult only).

It must be applied by qualified surgeons (orthopedists, neurosurgeons, cranio-maxillofacial surgeons, stomatologists and ENT surgeons) trained in bone grafting and fixation techniques and who have read these instructions for use.

The radio-opacity of BiologicGlass Injectable Putty makes it possible to discern bone substitute granules following their implantation. As the granules resorb, the radio-opacity of the bone defect approaches that of the surrounding bone.

BiologicGlass Injectable Putty biomaterial is a non-hardening paste, ready to use.

Contraindications

BiologicGlass Injectable Putty should not be used:

- In case of chronic or acute infection not treated with appropriate therapy.
- In patients who have suffered severe trauma with open external wounds near the defect, which are likely to become infected.
- In patients with a known allergy to bioactive glass or its constituents (Ca^{2+} , PO_4^{3-} , Na^+ and $\text{Si}(\text{OH})_4$), polyethylene glycol and / or glycerol.
- In patients with pre-existing conditions or disease that may interfere with proper tissue healing.
- In patients who have undergone or will undergo chemotherapy or radiation therapy at or near the site of implantation.
- In the irradiated bone (according to radiological criteria indicating osteonecrosis).
- To replace structures subjected to high mechanical stresses.
- During severe renal and hepatic infections.
- In conjunction with a treatment known to affect the skeleton.
- In case of unsutured meningeal breach in cranio-spinal surgery.
- In neonatology service.
- To date, we do not have any studies conducted in pregnant women or data related to use during breastfeeding.

As a safety measure, the implantation of BiologicGlass Injectable Putty is not recommended during the periods of pregnancy and lactation.

Composition

BiologicGlass Injectable Putty is an injectable paste composed of 45S5 bioactive glass granules (69%) mixed with a resorbable binder combining polyethylene glycol (12%) and glycerol (19%).

Mechanism of action / Performance / Benefits

During its implantation, BiologicGlass Injectable Putty is in contact with bone and biological fluids. Bioactive granules consist only of elements naturally present in bone tissue (Ca, P, Na, Si, O). After implantation, the binder is rapidly absorbed, allowing fluid circulation and cellular and vascular colonization of interstitial spaces between bioactive glass granules.

The resorption of bioactive glass will allow the formation on the surface of a layer of carbonate hydroxyapatite, which composition and structure are similar to the mineral phase of bone. This layer gives the granules their osteoconduction property and makes it possible to create a strong link between the granules and the living tissues.

The claimed clinical performance is the filling, reconstruction and / or fusion of bone defects allowing the regeneration of the bone.

The associated benefits are the absence/decrease of autologous bone sampling, a decrease in pain and an improvement in quality of life.

Instruction for use

- Check the expiry date. Do not use the product if it is exceeded.
- Do not use the device if the temperature indicator has turned or if the sterile packaging is damaged.
- Check each device before use, in order to detect any deterioration. If this is the case, do not use the implant.
- BiologicGlass Injectable Putty does not require mixing, humidification or other preparation prior to use. Humidification could dissolve the binder prematurely.
- Remove any soft and/or pathological tissues from the implantation site.
- Open the outer pouch (sterile barrier) and take out the inner pouch on the sterile field.
- Once the surgical site is prepared, open the internal pouch and remove the applicator.
- Unscrew the wing cap. Depending on the size of the defect and at the convenience of the surgeon, the tip can be unscrewed as well.
- Position the applicator at the level of the bone defect and press the plunger to gently fill the defect. Avoid direct contact between the delivery device and the patient's tissues. Use a sterile instrument if necessary. The device can also be placed directly in an intervertebral cage in spinal surgery.
- Fully fill the defect with BiologicGlass injectable Putty.
- It is also possible to deposit the paste in a sterile cup and then perform the implantation with a sterile instrument.
- The closure of the surgical site depends on the surgery performed and on the surgical site (membrane, stitches ...).
- This device is MR safe.

Warnings and precautions for use

In relation to the surgical procedure

- The general principles of asepsis and patient medication should be respected when using BiologicGlass Injectable Putty. BiologicGlass Injectable Putty does not substitute antibiotic therapy treatment during infection.
- The combination of any drug substance with BiologicGlass Injectable Putty during implantation is the responsibility of the surgeon.
- Handle BiologicGlass Injectable Putty with a surgical instrument to avoid piercing surgical gloves.
- It is advisable to trim up the recipient site before implantation.
- Avoid placing paste outside the bone defect. Remove it if necessary.
- If positioned outside the implantation site, moving or migrating, bioactive glass can cause wear of the joints and interfere with movement.
- Do not exert excessive pressure on the defect. Excessive pressure could cause an embolism of fat or paste in the bloodstream or cause the paste to be extruded beyond the implantation site, damaging the surrounding tissues.
- BiologicGlass Injectable Putty does not have sufficient mechanical strength to withstand a load before the bone tissue is formed. When used in load-bearing areas such as mandible fractures, standard internal or external stabilization techniques should be used to achieve rigid stabilization in all planes.
- It is necessary to follow the usual post-operative procedures of treatment and rehabilitation associated with bone grafts.

In relation to the medical device

- BiologicGlass Injectable Putty is a device that resorbs over time to make way for a regenerated bone. The binder is reabsorbed in a few days. Regarding the granules, no clinical studies currently available demonstrate complete resorption.
- BiologicGlass Injectable Putty is a non-hardening, ready-to-use paste.
- BiologicGlass Injectable Putty is a sterile single-use device and must not be re-sterilized or reused under any circumstances. Reuse can cause contamination and impaired performance of the bone substitute.

Adverse effects

No side effects directly related to the device have been reported to date. However, an unknown allergy to one of the constituents of the product could occur. Delayed union or failed fusion may also occur depending on the patient's metabolism.

Possible complications are general complications due to surgery or anesthesia: post-surgical symptoms (pain, redness, inflammation, edema, hematomas, seroma, swelling ...), postoperative infection, delay in consolidation, loss of fracture reduction, fusion failure, fracture, loss of bone graft, protuberance of the graft. These complications are the same as those that can occur with autologous bone grafts.

Any serious incident that may occur in connection with BiologicGlass Injectable Putty must be notified to NORAKER and the competent authority of the Member State in which the patient or surgeon is established.

Patient information

- The patient must be informed by the surgeon of the potential risks and adverse effects related to implantation and give his/her agreement to the proposed intervention.
- The surgeon must inform the patient receiving this device that the success of the implantation also depends on its behavior and good compliance with the post-operative hygiene instructions.
- The patient must report to his surgeon any event that may compromise the proper integration of the implant and undergo postoperative controls.
- After surgery, an implant card with its fascicle is filled out by the medical staff and given to the patient. He will have to keep it for life. Also, it is advisable to scan it when returning home.
- The SSCP (Summary of Safety and Clinical Performance) of the device is available on the manufacturer's website (www.noraker.com) (or on EUDAMED as soon as available).

Sterilization and packaging

Sterile device packed in sealed double-pouches and sterilized by gamma radiation (sterile barrier ensured by external pouch). Sterility is guaranteed until the expiry date if the sterile barrier has not been opened or damaged. Re-sterilization is prohibited.

Storage and disposal

Devices must be stored in their original unopened packaging, in a dry, clean place, protected from sunlight and at a recommended temperature between 15°C and 25°C.

The disposal of the device (packaging, syringe, remaining paste) must be carried out in accordance with local regulations and practices, at the risk of exposing users and patients to pathogens and contaminating the waste circuit.

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