# MedactaGlass

BIOACTIVE BONE SUBSTITUTE

EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute.....	5
FR – NOTICE D'UTILISATION - Substitut osseux bioactif .....	7
IT – ISTRUZIONI PER L'USO - Sostituto osseo bioattivo .....	9
DE – GEBRAUCHSANWEISUNG - Bioaktives Knochenersatzmaterial .....	12
ES - INSTRUCCIONES DE USO - Sustituto óseo bioactivo .....	14
PT - INSTRUÇÕES DE UTILIZAÇÃO - Substituto de osso bioativo .....	16
BG – ИНСТРУКЦИИ ЗА УПОТРЕБА - Биоактивен заместител на костна тъкан .....	18
CS – NÁVOD K POUŽITÍ - Bioaktivní kostní náhrada .....	21
EL – ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - Βιοενεργό υποκατάστατο οστών .....	23
HR – UPUTE ZA UPORABU - Bioaktivni koštani nadomjestak.....	25
NL – GEBRUIKSAANWIJZING - Bioactieve botvervanger.....	28
PL – INSTRUKCJA OBSŁUGI - Bioaktywny substytut kości.....	30
SL – NAVODILA ZA UPORABO - Bioaktivni nadomestek kosti .....	32
LV – LIETOŠANAS INSTRUKCIJA – Bioaktīvs kaulu aizstājējs .....	34

REF	VOL
XMG-IP1.0	1cc
XMG-IP2.5	2.5cc
XMG-IP5.0	5cc
XMG-IP10	10cc

Manufacturer :

**NORAKER**  
THE BIOGLASS® COMPANY


**NORAKER**  
60, Avenue Rockefeller  
69008 LYON, France  
Tel : +33 (0)4 78 93 30 92  
[contact@noraker.com](mailto:contact@noraker.com)  
[www.noraker.com](http://www.noraker.com)

Distributor :

**medacta**  
International 

**MEDACTA International**  
Strada Regina CH-6874  
Castel San Pietro, Switzerland  
Tel : +41 91 696 60 60  
[info@medacta.ch](mailto:info@medacta.ch)  
[www.medacta.ch](http://www.medacta.ch)

	FR	Conserver au sec
	EN	Keep dry
	IT	Tenere all'asciutto
	DE	Vor Feuchtigkeit schützen
	ES	Manténgase seco
	PT	Manter seco
	BG	Да се пази сухо
	CS	Uchovávejte v suchu
	EL	Να διατηρείται στεγνό
	HR	Čuvati na suhom
	NL	Droog houden
	PL	Trzymać w suchym miejscu
	SL	Hranite na suhem
	LV	Glabāt sausumā
	FR	Conserver à l'abri de la lumière du soleil
	EN	Store away from sunlight
	IT	Tenere lontano dalla luce del sole
	DE	Vor Sonnenlicht schützen
	ES	Mantener alejado de la luz solar
	PT	Manter afastado da luz solar
	BG	Да се пази от слънчева светлина
	CS	Uchovávejte mimo dosah slunečního záření
	EL	Να φυλάσσεται μακριά από το ηλιακό φως
	HR	Držati podalje od sunčeve svjetlosti
	NL	Uit de buurt van zonlicht houden
	PL	Trzymać z dala światła słonecznego
	SL	Hranite izven dosega sončne svetlobe
	LV	Glabāt prom no saules gaismas
	FR	Limite de température 15-25°C (stockage)
	EN	Temperature limit 15-25°C (storage)
	IT	Limite di temperatura 15-25 °C (conservazione)
	DE	Temperaturbegrenzung 15-25 °C (Lagerung)
	ES	Limite de temperatura 15-25 °C (almacenamiento)
	PT	Limite de temperatura 15-25°C (armazenamento)
	BG	Температурна граница 15-25°C (съхранение)
	CS	Teplotní omezení 15–25 °C (skladování)
	EL	Όριο θερμοκρασίας 15-25 °C (αποθήκευση)
	HR	Ograničenje temperature: 15-25 °C (skladištenje)
	NL	Temperatuurlimiet 15-25°C (opslag)
	PL	Limit temperatury 15–25°C (przechowywanie)
	SL	Temperaturna meja 15-25 °C (shranjevanje)
	LV	Temperatūras ierobežojums 15-25°C (uzglabāšana)
	AR	Ne pas réutiliser
	EN	Do not reuse
	IT	Non riutilizzare
	DE	Nicht erneut verwenden
	ES	No reutilizar
	PT	Não reutilizar
	BG	Да не се употребява повторно
	CS	Neroužívejte opakovaně
	EL	Μην επαναχρησιμοποιείτε
	HR	Nemojte ponovno upotrebljavati
	NL	Niet hergebruiken
	PL	Nie używać ponownie
	SL	Ponovna uporaba ni dovoljena
	LV	Nelietot atkārtoti
	FR	Stérile, Irradiation gamma
	EN	Sterile, Gamma irradiation
	IT	Sterile, irradiazione gamma
	DE	Sterilisiert durch Gamma-Bestrahlung
	ES	Estéril, radiación Gamma
	PT	Estéril, Irradiação gama
	BG	Стерилно, Гама облъчване
	CS	Sterilní, gama záření
	EL	Αποστειρωμένο με ακτινοβολία γάμμα
	HR	Sterilizirano gama-zračenjem
	NL	Steriel, gammastraling
	PL	Produkt sterylony, napromienianie gamma
	SL	Sterilno, gama obsevanje
	LV	Sterilizēts ar gamma starojumu
	FR	Ne pas restériliser
	EN	Do not resterilize
	IT	Non risterilizzare
	DE	Nicht erneut sterilisieren
	ES	No reesterilizar
	PT	Não reesterilizar
	BG	Да не се стерилизира повторно
	CS	Nesterilizujte opakovaně
	EL	Μην επαναποστεριώνετε
	HR	Nemojte ponovno sterilizirati
	NL	Niet opnieuw steriliseren
	PL	Nie sterylizować ponownie
	SL	Ne sterilizirajte ponovno
	LV	Nesterilizēt atkārtoti
FR	Ne pas utiliser si l'emballage est endommagé et consulter la notice d'utilisation	

	EN	Do not use if package is damaged and consult instruction for use
	IT	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	DE	Nicht verwenden, wenn die Verpackung beschädigt ist und Gebrauchsanweisung beachten
	ES	No utilizar si el envase está dañado. Consultense las instrucciones de uso.
	PT	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	BG	Да не се употребява, ако опаковката е повредена и да се следват инструкциите за употреба
	CS	Pokud je obal poškozen, nepoužívejte a přečtěte si návod k použití
	EL	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευθείτε τις οδηγίες χρήσης
	HR	Nemojte upotrebljavati ako je pakiranje oštećeno; pogledajte upute za uporabu
	NL	Niet gebruiken als de verpakking beschadigd is. Raadpleeg de gebruiksaanwijzing
	PL	Nie używać w razie uszkodzenia opakowania i zapoznać się z instrukcją obsługi
	SL	Ne uporabljajte, če je ovojnina poškodovana. Oglejte si navodila za uporabo
	LV	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju
	FR	Attention Notice d'utilisation
	EN	Instruction for Use
	IT	Attenzione Istruzioni per l'uso
	DE	Vorsicht
	ES	Precaución, Instrucciones de uso
	PT	Instruções de cuidado de utilização
	BG	Внимание Инструкции за употреба
	CS	Upozornění: Návod k použití
	EL	Προσοχή, συμβουλευτείτε τις οδηγίες χρήσης
	HR	Oprez
	NL	Opgelet gebruiksaanwijzing
	PL	Przeostroga – instrukcja obsługi
	SL	Navodila za uporabo
	LV	Lietošanas instrukcija
	FR	Numéro de lot
	EN	Batch number
	IT	Numero di lotto
	DE	Chargenbezeichnung
	ES	Código de lote
	PT	Numero do lote
	BG	Номер на партида
	CS	Číslo dávky
	EL	Αριθμός партиδας
	HR	Broj serije
	NL	Partijnummer
	PL	Numer partii
	SL	Številka serije
	LV	Partijas numurs
	FR	Date limite d'utilisation
	EN	Use by date
	IT	Data di scadenza
	DE	Verwendbar bis
	ES	Fecha de caducidad
	PT	Data de validade
	BG	Срок на годност
	CS	Použitelné do
	EL	Ημερομηνία λήξης
	HR	Rok uporabe
	NL	Houdbaarheidsdatum
	PL	Termin przydatności
	SL	Datum uporabe
	LV	Izlietot līdz (datums)
	FR	Référence
	EN	Reference
	IT	Riferimento
	DE	Bestellnummer/Artikelnummer
	ES	Referencia
	PT	Referência
	BG	Референция
	CS	Reference
	EL	Αναφορά
	HR	Ref. oznaka
	NL	Referentie
	PL	Odnosnik
	SL	Reference
	LV	Atsauce
		FR
EN		Device marked CE by the GMED notified body. Placed on the market in 2017
IT		Dispositivo marcato CE dall'organismo notificato GMED. Immeso sul mercato nel 2017
DE		Produkt mit CE-Kennzeichnung durch die benannte Stelle GMED. Markteinführung 2017

	ES	Producto con marcado CE del organismo notificado GMED. Comercializado desde 2017
	PT	Dispositivo com a marca CE emitida pelo organismo notificado pela GMED. Lançado no mercado em 2017
	BG	Изделие с ЕС-знак от GMED нотифициран орган. Пуснато на пазара от 2017 г.
	CS	Zdravotnícký prostředek s označením EC notifikovaným orgánem GMED. Uvedeno na trh v roce 2017
	EL	Προϊόν με σήμανση EC από τον κοινοποιημένο οργανισμό GMED. Διατέθηκε στην αγορά το 2017
	HR	Prijavljeno tijelo GMED označilo je proizvod oznakom EC. Stavljeno na tržište 2017.
	NL	Apparaat EG-gemarkeerd door de GMED-aangemelde instantie. Op de markt gebracht in 2017
	PL	Wyrób oznaczony symbolem WE przez jednostkę notyfikowaną GMED. Rok wprowadzenia na rynek: 2017
	SL	Pripomoček ES – označen s strani priglašenega organa GMED. Na trg je bil dan leta 2017
	LV	Ierīcei ir CE marķējums, ko veikusi tās pilnvarotā iestāde GMED. Laista tirgū 2017. gadā
	FR	Fabricant
	EN	Manufacturer
	IT	Produttore
	DE	Hersteller
	ES	Fabricante
	PT	Fabricante
	BG	Производител
	CS	Výrobce
	EL	Κατασκευαστής
	HR	Proizvođač
	NL	Fabrikant
	PL	Producent
	SL	Proizvajalec
LV	Ražotājs	
		NORAKER 60, Avenue Rockefeller 69008 Lyon, FRANCE Tel : +33 (0)4 78 93 30 92 Fax : +33 (0)4 72 35 94 37
	FR	Date de fabrication
	EN	Manufacturing date
	IT	Data di produzione
	DE	Herstellungsdatum
	ES	Fecha de fabricación
	PT	Data de fabrico
	BG	Дата на производство
	CS	Datum výroby
	EL	Ημερομηνία κατασκευής
	HR	Datum proizvodnje
	NL	Productiedatum
	PL	Data produkcji
	SL	Datum proizvodnje
LV	Ražošanas datums	
	FR	Dispositif Medical
	EN	Medical Device
	IT	Dispositivo medico
	DE	Medizinprodukt
	ES	Producto sanitario
	PT	Dispositivo Médico
	BG	Медицинско изделие
	CS	Zdravotnícký prostředek
	EL	Ιατροτεχνολογικό προϊόν
	HR	Medicinski proizvod
	NL	Medisch apparaat
	PL	Wyrób medyczny
	SL	Medicinski pripomoček
LV	Medicīnas ierīce	
	FR	Consulter les instructions d'utilisation
	EN	Consult instructions for use
	IT	Consultare le istruzioni per l'uso
	DE	Siehe Gebrauchsanweisung
	ES	Consultense las instrucciones de uso
	PT	Consultar as instruções de utilização
	BG	Следвайте инструкциите за употреба
	CS	Přečtěte si návod k použití
	EL	Συμβουλευθείτε τις οδηγίες χρήσης
	HR	Obvezno pročitati upute
	NL	Raadpleeg de gebruiksaanwijzing
	PL	Zapoznać się z instrukcją obsługi
	SL	Ravnajte se po navodilih za uporabo
LV	Skatīt lietošanas instrukciju	
	FR	Volume (cc / cm <sup>3</sup> )
	EN	Volume (cc / cm <sup>3</sup> )
	IT	Volume (cc / cm <sup>3</sup> )
	DE	Volumen (cm <sup>3</sup> )
	ES	Volumen (cm <sup>3</sup> )
	PT	Volume (cc / cm <sup>3</sup> )
	BG	Обем (cc / cm3)
	CS	Objem (cc/cm3)
	EL	Όγκος (cc/cm3)
	HR	Volumen (cc / cm3)
	NL	Volume (cc/cm3)
	PL	Pojemność (cc / cm3)

	SL	Prostornina (cc/cm3)
	LV	Tilpums (cc / cm3)
	FR	Système de barrière stérile unique avec emballage de protection intérieure
	EN	Single sterile barrier system with protective packaging inside
	IT	Sistema a barriera sterile singola con confezione protettiva interna
	DE	Einfach-Steril-Barrieresystem mit inliegender Schutzverpackung
	ES	Sistema protector estéril único con envase de protección interno
	PT	Sistema de barreira estéril única com embalagem protetor no interior
	BG	Система с единична стерилна бариера с вътрешна защитна опаковка
	CS	Jediný sterilní bariérový systém s ochranným obalem uvnitř
	EL	Ενιαίο σύστημα στείρου φραγμού με εσωτερική προστατευτική συσκευασία
	HR	Sustav s jednostrukom sterilnom barijerom sa zaštitnim unutarnjim pakiranjem
	NL	Enkel steriel barrièresysteem met beschermende verpakking aan de binnenkant
	PL	System pojedynczej bariery sterylnej z wewnętrznym opakowaniem ochronnym
	SL	Sistem sterilne pregrade z zaščitno embalažo v notranjosti
	LV	Viena sterila barjersistēma ar aizsargiepakojumu iekšpusē
	FR	Identifiant unique de dispositif
	EN	Unique Device Identifier
	IT	Identificativo univoco del dispositivo
	DE	Einmalige Produktkennung
	ES	Identificador único del dispositivo (UDI)
	PT	Identificador Único do Dispositivo
	BG	Уникален идентификатор на изделието
	CS	Jedinečný identifikátor zařízení
	EL	Μοναδικό αναγνωριστικό προϊόντος
	HR	Jedinstvena identifikacija proizvoda
	NL	Unieke apparaat-ID
	PL	Unikalny identyfikator urządzenia
	SL	Enolični identifikator naprave
	LV	Ierīces unikālais identifikators
<b>2023/08</b>	FR	Mise à jour du document
	EN	Last update
	IT	Ultimo aggiornamento
	DE	Letzte Überarbeitung
	ES	Última actualización
	PT	Última atualização
	BG	Последна промяна
	CS	Poslední aktualizace
	EL	Τελευταία ενημέρωση
	HR	Zadnji put revidirano
	NL	Laatste update
	PL	Ostatnia aktualizacja
	SL	Zadnja posodobitev
	LV	Pēdējais atjauninājums

## EN – INSTRUCTIONS FOR USE – Bioactive Bone substitute

### Description and indications

MectaGlass Injectable Putty is a synthetic, resorbable, biocompatible and bioactive substitute device (bioactive glass 45S5), intended for the filling, reconstruction and / or fusion of bone defects or gaps of the skeletal system, in orthopedic surgery, spine, cranio-maxillofacial surgery and ENT.

It is indicated in case of loss or lack of bone substance for bone defects of traumatic, pathological or surgical origin when autologous solutions are not applicable or sufficient in orthopedics, neurosurgery, cranio maxillo facial and otorhinolaryngology surgery in children and adult population:

- Fusion or reconstruction of deformities and degenerative diseases in spine
- Fusion or reconstruction of deformities and degenerative bone pathologies in orthopedic
- Filling and reconstruction of bone defects due to resection of tumors, cyst or infection and in case of prosthetic revision
- Filling after surgical bone defect (donor sites after removal of autograft, trepanation,...)
- Filling after removal of cholesteatoma
- Filling and reconstruction due to maxilla and periodontium pathologies (in adult only).

It must be applied by qualified surgeons (orthopedists, neurosurgeons, cranio-maxillofacial surgeons, stomatologists and ENT surgeons) trained in bone grafting and fixation techniques and who have read these instructions for use.

The radio-opacity of MectaGlass Injectable Putty makes it possible to discern bone substitute granules following their implantation. As the granules resorb, the radio-opacity of the bone defect approaches that of the surrounding bone.

MectaGlass Injectable Putty biomaterial is a non-hardening paste, ready to use.

### Contraindications

MectaGlass Injectable Putty should not be used:

- In case of chronic or acute infection not treated with appropriate therapy.
- In patients who have suffered severe trauma with open external wounds near the defect, which are likely to become infected.

- In patients with a known allergy to bioactive glass or its constituents (Ca<sup>2+</sup>, PO<sub>4</sub><sup>3-</sup>, Na<sup>+</sup> and Si (OH)<sub>4</sub>), polyethylene glycol and / or glycerol.
- In patients with pre-existing conditions or disease that may interfere with proper tissue healing.
- In patients who have undergone or will undergo chemotherapy or radiation therapy at or near the site of implantation.
- In the irradiated bone (according to radiological criteria indicating osteonecrosis).
- To replace structures subjected to high mechanical stresses.
- During severe renal and hepatic infections.
- In conjunction with a treatment known to affect the skeleton.
- In case of unsutured meningeal breach in cranio-spinal surgery.
- In neonatology service.
- To date, we do not have any studies conducted in pregnant women or data related to use during breastfeeding.

As a safety measure, the implantation of MectaGlass Injectable Putty is not recommended during the periods of pregnancy and lactation.

#### Composition

MectaGlass Injectable Putty is an injectable paste composed of 45S5 bioactive glass granules (69%) mixed with a resorbable binder combining polyethylene glycol (12%) and glycerol (19%).

#### **Mechanism of action / Performance / Benefits**

During its implantation, MectaGlass Injectable Putty is in contact with bone and biological fluids. Bioactive granules consist only of elements naturally present in bone tissue (Ca, P, Na, Si, O). After implantation, the binder is rapidly absorbed, allowing fluid circulation and cellular and vascular colonization of interstitial spaces between bioactive glass granules.

The resorption of bioactive glass will allow the formation on the surface of a layer of carbonate hydroxyapatite, which composition and structure are similar to the mineral phase of bone. This layer gives the granules their osteoconduction property and makes it possible to create a strong link between the granules and the living tissues.

The claimed clinical performance is the filling, reconstruction and / or fusion of bone defects allowing the regeneration of the bone.

The associated benefits are the absence/decrease of autologous bone sampling, a decrease in pain and an improvement in quality of life.

#### **Instruction for use**

- Check the expiry date. Do not use the product if it is exceeded.
- Do not use the device if the temperature indicator has turned or if the sterile packaging is damaged.
- Check each device before use, in order to detect any deterioration. If this is the case, do not use the implant.
- MectaGlass Injectable Putty does not require mixing, humification or other preparation prior to use. Humidification could dissolve the binder prematurely.
- Remove any soft and/or pathological tissues from the implantation site.
- Open the outer pouch (sterile barrier) and take out the inner pouch on the sterile field.
- Once the surgical site is prepared, open the internal pouch and remove the applicator.
- Unscrew the wing cap. Depending on the size of the defect and at the convenience of the surgeon, the tip can be unscrewed as well.
- Position the applicator at the level of the bone defect and press the plunger to gently fill the defect. Avoid direct contact between the delivery device and the patient's tissues. Use a sterile instrument if necessary. The device can also be placed directly in an intervertebral cage in spinal surgery.
- Fully fill the defect with MectaGlass injectable Putty.
- It is also possible to deposit the paste in a sterile cup and then perform the implantation with a sterile instrument.
- The closure of the surgical site depends on the surgery performed and on the surgical site (membrane, stitches ...).
- This device is MR safe.

#### **Warnings and precautions for use**

##### ***In relation to the surgical procedure***

- The general principles of asepsis and patient medication should be respected when using MectaGlass Injectable Putty. MectaGlass Injectable Putty does not substitute antibiotic therapy treatment during infection.
- The combination of any drug substance with MectaGlass Injectable Putty during implantation is the responsibility of the surgeon.
- Handle MectaGlass Injectable Putty with a surgical instrument to avoid piercing surgical gloves.
- It is advisable to trim up the recipient site before implantation.
- Avoid placing paste outside the bone defect. Remove it if necessary.
- If positioned outside the implantation site, moving or migrating, bioactive glass can cause wear of the joints and interfere with movement.
- Do not exert excessive pressure on the defect. Excessive pressure could cause an embolism of fat or paste in the bloodstream or cause the paste to be extruded beyond the implantation site, damaging the surrounding tissues.
- MectaGlass Injectable Putty does not have sufficient mechanical strength to withstand a load before the bone tissue is formed. When used in load-bearing areas such as mandible fractures, standard internal or external stabilization techniques should be used to achieve rigid stabilization in all planes.

- It is necessary to follow the usual post-operative procedures of treatment and rehabilitation associated with bone grafts.

#### **In relation to the medical device**

- MectaGlass Injectable Putty is a device that resorbs over time to make way for a regenerated bone. The binder is reabsorbed in a few days. Regarding the granules, no clinical studies currently available demonstrate complete resorption.
- MectaGlass Injectable Putty is a non-hardening, ready-to-use paste.
- MectaGlass Injectable Putty is a sterile single-use device and must not be re-sterilized or reused under any circumstances. Reuse can cause contamination and impaired performance of the bone substitute.

#### **Adverse effects**

No side effects directly related to the device have been reported to date. However, an unknown allergy to one of the constituents of the product could occur. Delayed union or failed fusion may also occur depending on the patient's metabolism.

Possible complications are general complications due to surgery or anesthesia: post-surgical symptoms (pain, redness, inflammation, edema, hematomas, seroma, swelling ...), postoperative infection, delay in consolidation, loss of fracture reduction, fusion failure, fracture, loss of bone graft, protuberance of the graft. These complications are the same as those that can occur with autologous bone grafts.

Any serious incident that may occur in connection with MectaGlass Injectable Putty must be notified to NORAKER and the competent authority of the Member State in which the patient or surgeon is established.

#### **Patient information**

- The patient must be informed by the surgeon of the potential risks and adverse effects related to implantation and give his/her agreement to the proposed intervention.
- The surgeon must inform the patient receiving this device that the success of the implantation also depends on its behavior and good compliance with the post-operative hygiene instructions.
- The patient must report to his surgeon any event that may compromise the proper integration of the implant and undergo postoperative controls.
- After surgery, an implant card with its fascicle is filled out by the medical staff and given to the patient. He will have to keep it for life. Also, it is advisable to scan it when returning home.
- The SSCP (Summary of Safety and Clinical Performance) of the device is available on the manufacturer's website ([www.noraker.com](http://www.noraker.com)) (or on EUDAMED as soon as available).

#### **Sterilization and packaging**

Sterile device packed in sealed double-pouches and sterilized by gamma radiation (sterile barrier ensured by external pouch). Sterility is guaranteed until the expiry date if the sterile barrier has not been opened or damaged. Re-sterilization is prohibited.

#### **Storage and disposal**

Devices must be stored in their original unopened packaging, in a dry, clean place, protected from sunlight and at a recommended temperature between 15°C and 25°C.

The disposal of the device (packaging, syringe, remaining paste) must be carried out in accordance with local regulations and practices, at the risk of exposing users and patients to pathogens and contaminating the waste circuit.

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