Activloss®

Bioactive Bone Substitute

GRANULES

- EN - Caution Instruction for Use
- AR - تديه تعليمات الاستخدام
- BG - Внимание Инструкция за употреба
- DE - Vorsicht
- ES - Precaución, instrucciones de uso
- FR - Attention Notice d'utilisation
- HR - Oprez
- HU - Vigyázat – Használati utasítás
- IT - Attenzione Istruzioni per l'uso
- LV - Uzmanību Lietošanas instrukcija
- PT - Instruções de cuidado de utilização
- SL - Previdnost Navodila za uporabo
- SR - Oprez Упутство за употребу

NORAKER
60, Avenue Rockefeller
69008 LYON - FRANCE
Tel : +33 (0)4 78 93 30 92

- EN : Manufacturer
- AR : الشركة المصنعة
- BG : Производител
- DE : Hersteller
- ES : Fabricante
- FR : Fabricant
- HR : Proizvođač
- HU : Valmistaja
- IT : Produttore
- LV : Ražotājs
- PT : Fabricante
- SL : Proizvajalec
- SR : Proizvođač

REF	 (mm)	VOL (cc/cm3)
EN : References	EN : Granule size	EN : Volume
AR : المرجع	AR : قيس الحبيبات	AR : الحجم (ستيمتر مكعب)
BG : Референция	BG : Размер на гранулите	BG : Обем
DE : Bestellnummern	DE : Korngröße	DE : Volumen
ES : Referencia	ES : Granulometría	ES : Volumen
FR : Références	FR : Granulométrie	FR : Volume
HR : Ref. oznaka	HR : Veličina granula	HR : Volumen
HU : Hivatkozásk	HU : Szecseméret	HU : Térfogat
IT : Riferimenti	IT : Granulometria	IT : Volume
LV : Atsaucē	LV : Granulu izmērs	LV : Tilpums
PT : Referências	PT : Granulometria	PT : Volume
SL : Referenca	SL : Velikost zrnč	SL : Količina
SR : Референтна ознаке	SR : Величина гранула	SR : Запремина
ACT-GM0.5		0,5cc
ACT-GM1	0.5-1mm / 1 mm	1cc
ACT-GM2		2cc
ACT-GL1	1-3mm / 3 mm	1cc

EN Consult instructions for use	اطلع على تعليمات الاستخدام
AR	اقرأ التعليمات الاستخدام
BG	Следвайте инструкциите за употреба
DE	Gebräuchsanweisung beachten
ES	Consulte las instrucciones de uso
FR	Consulter les instructions d'utilisation
HR	Obvezno pročitati upute
HU	Olvassa el a használati utasítást
IT	Consultare le istruzioni per l'uso
LV	Skatīt lietošanas instrukciju
PT	Consultar as instruções de utilização
SL	Ravnajte se po navodilih za uporabo
SR	Pogledati uputstva za upotrebu

EN Do not use if package is damaged and consult instruction for use	لا تستخدم المنتج في حالة تلف العبوة واطلع على تعليمات الاستخدام
AR	لا تستخدم المنتج في حالة تلف العبوة واطلع على تعليمات الاستخدام
BG	Да не се употребява, ако опаковката е повредена и да се следват инструкциите за употреба
DE	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten
ES	No utilizar si el envase está dañado y consúltense las instrucciones de uso.
FR	Ne pas utiliser si l'emballage est endommagé et consulter la notice d'utilisation
HR	Nemojte upotrebljavati ako je pakiranje oštećeno; pogledajte upute za upotrebu
HU	Ne használja, ha a csomagolás sérült, és olvassa el a használati utasítást
IT	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
LV	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju
PT	Não utilizar se a embalagem estiver danificada; consultar as instruções de utilização
SL	Ne uporabljajte, če je ovojnina poškodovana. Oglejte si navodila za uporabo
SR	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu

EN Device EC-marked by the GMED notified body. Placed on the market in 2008	الجهاز معتمد طرُح في GMED معتمد بعلامة المنظمة الأوروبية من قبل جهة السوق في عام 2008
BG	Изделие с ЕС-знак от GMED нотифициран орган. Пуснато на пазара от 2008 г.
DE	Produkt mit CE-Kennzeichnung durch die benannte Stelle GMED. Markteinführung 2008
ES	Producto con marcado CE del organismo notificado GMED. Comercializado desde 2008
FR	Dispositif marqué CE par l'organisme notifié GMED. Mis sur le marché en 2008
HR	Prijavljeno tijelo GMED označilo je proizvod oznakom EC. Stavljeno na tržište 2008
HU	A GMED bejelentett szervezet által EC-jelöléssel ellátott eszköz. Forgalomba hozva 2008-ban

EN	Device EC-marked by the GMED notified body. Placed on the market in 2008
AR	المعتمد طرُح في GMED معتمد بعلامة المنظمة الأوروبية من قبل جهة السوق في عام 2008
BG	Изделие с ЕС-знак от GMED нотифициран орган. Пуснато на пазара от 2008 г.
DE	Produkt mit CE-Kennzeichnung durch die benannte Stelle GMED. Markteinführung 2008
ES	Producto con marcado CE del organismo notificado GMED. Comercializado desde 2008
FR	Dispositif marqué CE par l'organisme notifié GMED. Mis sur le marché en 2008
HR	Prijavljeno tijelo GMED označilo je proizvod oznakom EC. Stavljeno na tržište 2008
HU	A GMED bejelentett szervezet által EC-jelöléssel ellátott eszköz. Forgalomba hozva 2008-ban

IT	Dispositivo marcato CE dall'organismo notificato GMED. Imnesso sul mercato nel 2008
LV	Ierīcei ir EK marķējums, ko piešķirusi GMED pilnvarotā iestāde. Laista tirgū 2008. gadā.
PT	Dispositivo com a marca CE emitida pelo organismo notificado pela GMED. Lançado no mercado em 2008
SL	Priporoček ES in označen s strani prijavljenega organa GMED. Na trg je bil dan leta 2008
SR	Средство означено ознаком ЕС од стране GMED нотификованог тела. Стављено на тржиште 2008. год.

EN	Last update
AR	تاريخ آخر تحديث
BG	Последна промяна
DE	Letzte Überarbeitung
ES	Última actualización
FR	Mise à jour du document
HR	Zadnji put revidirano
HU	Legutóbbi frissítés
IT	Ultimo aggiornamento
LV	Pēdējais atjauninājums
PT	Última atualização
SL	Zadnja posodobitev
SR	Последње ажурирање

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EN – INSTRUCTION FOR USE

Bioactive bone substitute

Description and Indications
Activioss is a synthetic and biocompatible bone substitute device (bioactive glass 45S5), intended for the filling or reconstruction of bone defects or bone preservation in dental and maxillofacial systems surgery.
Activioss is indicated in case of loss or lack of bone substance for bone defects of surgical, traumatic or pathological origin when autologous solutions are not applicable or sufficient in adult population:
- Limited elevation and filling before implantation
- Periodontal defects filling and reconstruction
- Extraction sites, (alveolar ridge preservation, implant preparation, reconstruction or augmentation)
- Cyst cavities filling

It must be used by qualified surgeons (dentist surgeons, implantologists and stomatologists) trained in bone grafting and fixation techniques who have read these instructions for use.

Contraindications
Activioss should not be used:
- In case of chronic or acute infection not treated with appropriate therapy
- In patients who have suffered serious trauma with external wounds open near the defect, which could become infected.
- In patients with known allergy to bioactive glass or its constituents (Ca²⁺, PO³⁻, Na⁺ and Si(OH)₄).
- In patients with pre-existing conditions or disease that may interfere with the good healing of tissues
- In the irradiated bone (according to radiological criteria indicating osteonecrosis) and in special clinical situations (tumor, ongoing chemotherapy and radiation -therapy, immunodeficiency ...)
- During severe renal and hepatic infections.
In conjunction with a treatment known to affect the skeleton.
To date, we do not have any studies conducted in pregnant women or data related to use during breastfeeding. As a safety measure, the implantation of Activioss is not recommended during the periods of pregnancy and lactation.

Composition of Activioss

Activioss (45S5) consist only of elements found naturally in bone tissue: 24.5% calcium oxide (CaO), 6% phosphate pentoxide (P₂O₅), 24.5% sodium oxide (Na₂O) and 45% silicon dioxide (SiO₂).

Mechanism of action / Performance / Benefits

During its implantation, Activioss is in contact with bone and biological fluids. The release of ions during resorption will allow the surface formation of a layer of carbonate hydroxyapatite, whose composition and structure are similar to the mineral phase of bone. This layer gives Activioss its osteoconduction property and creates a strong link between the granules and living tissue.
The claimed clinical performance is the filling or reconstruction of bone defects allowing the regeneration of the bone.
The main associated benefit is the absence of morbidity (pain and complication) associated with harvesting graft.
Additional benefits with quality-of-life improvement are also identified:

-Maintain/survival implant (if the implant is used)
-Limited Recurrent/ Residual disease (for cyst resection)
-Optimize the esthetical and functional outcome (e.g, comfort, speech, chewing), visit decrease, easy to use for surgeons

All these benefits (main and additional) influence the quality-of-life improvement of Activioss patients.

Instructions for use
-Check the expiry date. Do not use the product if it is exceeded.
-Do not use if the sterile packaging (external blister) is damaged, opened prematurely or exposed to environmental conditions other than those specified.
-Check each device before use to detect any deterioration. If deterioration is present, do not use the implant.
-Eliminate all soft and/or pathological tissue from the implantation site.
-Open the outer blister (sterile barrier) and remove the inner blister on the sterile field.

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آلية العمل / الأداء / الفوائد
-The combination of any drug substance with Activioss during implantation is the responsibility of the surgeon.
-Manipulate Activioss with a surgical instrument to avoid piercing surgical gloves.
-It is advisable to revise the recipient site before implantation.
-Completely fill the defect with Activioss. It is possible to perform the application of Activioss if the defect has sufficient bone wall.
-Avoid placing granules outside of the bone defect. Remove them if necessary.
-Avoid direct contact of Activioss with the skin.
-If it moves/migrates, the bioactive glass can cause wear of the joints and interfere with movement. Prevention of movement and granule migration is essential for proper bone formation.
-Do not apply excessive pressure to the defect. Excessive pressure may cause embolization of fat in the bloodstream.
-Activioss maintains its volume that is to say it does not shrink or expand.
-Activioss does not have sufficient mechanical strength to withstand load bearing before hard tissue is formed. When used in load bearing areas such as mandible fractures, standard internal or external stabilization techniques should be followed to achieve rigid stabilization in all planes.
-Activioss should not be used in applications with immediate load bearing. Placing a dental implant under load can be done 5 to 6 months after filling. Immediate loading is the responsibility of the surgeon.
-It is necessary to follow the usual post-operative treatment and rehabilitation procedures associated with bone grafts.
-The closure of the operative site depends on the surgery performed and the surgical site (membrane, sutures, etc.). An adequate closure of the graft site is mandatory (e.g., collagen or synthetic membrane (resorbable or non-resorbable), mucosal-periosteum flap).

-The use of a membrane is recommended in cases where the defect is significant or if there is limited bone retention. Primary closure of defects, preferably without tension, is strongly recommended.

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Regarding the medical device
-Activioss is a device that resorbs over time to make way for regenerated bone.
-There is currently no clinical study available that demonstrates complete resorption of the granules.
-This device does not harden like cement.

-Activioss is a sterile disposable device and must never be re-sterilized or reused. Reuse may cause contamination and impairment of bone substitute performance.

Adverse effects

No side effect directly linked to the device has been reported to date. However, an unknown allergy to one of the constituents of the product may be possible. Possible complications or adverse events are not more severe than those expected of similar products if the instructions are followed correctly by a qualified surgeon familiar with bone grafting techniques: post-surgical symptoms (pain, discomfort, abscess, redness, inflammation, oedema, hematomas, seroma, swelling, bleeding, ...), postoperative infection, recurrence/residual disease, wound dehiscence, wound leakage, delay in consolidation, loss of bone graft, membrane exposure (if applicable), protuberance of the graft, tooth sensitivity, gingival recession, subsidence of the flap, abscess formation, resorption or ankylosis of the treated root, etc.

These complications are the same as those that can occur with the autologous bone grafting.

Any serious incident which may occur in connection with Activioss must be notified to NORAKER and to the competent authority of the Member State in which the patient or the surgeon is established.

Patient information

-The patient must be informed by the surgeon of the potential risks and adverse effects of bone substitution and agrees to the proposed procedure.
-The surgeon should inform the patient who is the recipient of this device that the success of bone substitution depends on the behavior and good compliance with post-operative hygiene instructions.

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